



The Republic of Uganda
Ministry of Health

Health Infrastructure Department

Operation Manual for
Regional Medical Equipment Maintenance
Workshops and
Medical Equipment Maintenance Guidelines

Volume I

A guide for Engineers and Technicians working on
Medical Equipment in Health Facilities

August 2020

Table of Contents

Foreword	5
Acronyms	7
Chapter 1: Introduction	9
1.1 Background	
1.2 Overall Objectives	
1.3 Specific Objectives	
1.4 Stakeholders in Medical Equipment Management	
1.5 Definition of Medical Equipment Maintenance	
Chapter 2: Organization Structure for Medical Equipment Maintenance	14
2.1 NACME	
2.2 Health Infrastructure Department (HID)	
2.3 Regional Medical Equipment Maintenance Workshops (RWs)	
2.4 National Referral Hospitals (NRHs)	
2.5 Regional Referral Hospitals (RRHs)	
2.6 District Health Services	
Chapter 3: Function, Operation and Role of RWs	20
3.1 Central Medical Equipment Maintenance Workshop	
3.2 Regional Medical Equipment Maintenance Workshops	
3.3 RW Maintenance Work Flow Chart	

Chapter 4: Planning, Budgeting and Management of Medical Equipment Maintenance 25

- 4.1 Medical Equipment Inventory Collection and Update
 - 4.1.1 Inventory taking Rules and Procedure for Inventory Update
 - 4.1.2 Basic Medical Equipment Inventory Data Analysis and Use
- 4.2 Work Plan and Budgeting for Medical Equipment Maintenance
 - 4.2.1 Process of Preparation of Annual Work Plans and Budgets
 - 4.2.2 Preparation of RW Annual and Quarterly Budget
 - 4.2.3 Preparation of Quarterly Maintenance Schedules
 - 4.2.4 Budgeting for Emergency Works
- 4.3 Medical Equipment Maintenance Planning and Management
 - 4.3.1 Preparation of Quarterly Reports
 - 4.3.2 RWs Performance Review Meetings
 - 4.3.3 Regional Workshop Management Committee (RWMC)
 - 4.3.4 Internal Workshop Staff Meetings at Workplace

Chapter 5: Procurement and Disposal of Goods and Services by RWs 38

- 5.1 Management of Procurement by RWs
 - 5.1.1 Determination of Spare Parts Needs
 - 5.1.2 Procedure for Purchasing Spare Parts
 - 5.1.3 Stores Management by RWs
- 5.2 Disposal of Items

Chapter 6:	Support Supervision, Monitoring and Evaluation of RWs	43
6.1	Support Supervision for RWs	
Chapter 7:	Guidelines for Medical Equipment Maintenance	45
7.1	Introduction	
7.2	Planned Preventive Maintenance (PPM)	
7.3	Corrective Maintenance Services	
7.4	Maintenance Policy Guideline	
7.5	Levels of Maintenance	
7.6	Medical Equipment that can be maintained by RWs or by Outsourcing the Services	
7.7	Setting up a Maintenance System	
7.7.1	Routine Maintenance Procedures	
7.7.2	Engineering Personnel	
7.7.3	Reminder System	
7.7.4	Surveillance	
7.7.5	Standard Maintenance Work Formats	
7.7.6	Special Test Equipment	
7.7.7	Technical Library	
Chapter 8:	Periodic Maintenance Checklist for Common Equipment (31 types of Equipment)	56

- Annex 1: Standard Floor Plan of a Regional Workshop
- Annex 2: Recommended RW Tools, Test Equipment and Furniture
- Annex 3: Equipment Condition A–F
- Annex 4: Medical Equipment Inventory Form
- Annex 5: Annual Work Plan Form
- Annex 6(a): 5S Activity’s Instruction for RWs
- Annex 6(b): CQI Action Plan Form
- Annex 7: Quarterly Report Form
- Annex 8: Job Card Form
- Annex 9: Annual Assessment Sheet for RWs
- Annex 10: Maintenance Requisition Form
- Annex 11: Maintenance Service Sticker Form

Foreword

The mission of the Ministry of Health is “to provide the highest possible level of health services to all people in Uganda through delivery of promotive, preventive, curative, palliative and rehabilitative health services at all levels”. In pursuit of the above mission, government ensures that healthcare services are delivered in an efficient, safe and cost-effective manner.

In this regard, a substantial proportion of the Health Sector budget is spent on acquisition of Medical Equipment. To ensure that available Medical Equipment safely works for a long time, it must be used and maintained properly to guarantee its reliability and the quality of Healthcare delivered to patients.

Appropriate preventive and corrective maintenance are key to achieving safe and cost-effective management of medical equipment. It is therefore, important that measures are taken to ensure medical equipment is well maintained and cared for by the Healthcare workers in order to maximize the investment made in its acquisition. However, for effective maintenance to be carried out by the maintenance teams, adequate funds should be budgeted and allocated.

This manual helps to streamline the operations of the Regional Medical Equipment Maintenance Workshops (RWs). Additionally, the Maintenance guidelines will facilitate first line Maintenance to be undertaken by both the users of the equipment and Engineers/Technicians in a well-structured and safe manner. Additionally, it is meant to be a reference book for the RWs to plan, manage and execute proper maintenance on commonly used essential Medical Equipment maintenance using the guidelines herein.

The SOPs for carrying out PPM (Volume II), will enable the Engineers and Technicians to maintain high level of work consistency and safety; resulting into equipment reliability and uptime.

I am pleased that the second edition of the manual has been updated to enable Engineers and Technicians working in health facilities and RWs to cope with the current challenges of managing medical equipment including electrical safety testing and calibration. I also recommend the use of this manual to carry out User Training for health workers.

Special recognition goes to Mr. Naoki Mimuro, Eng. Sitra Mulepo C.S and Eng. Owen Muhimbise for their dedication and Technical input to produce this manual.

Lastly, I would like to extend my appreciation to the staff of the Health Infrastructure Department and Regional Medical Equipment Maintenance Workshops, Members of National Advisory Committee on Medical Equipment (NACME) and Health Infrastructure Technical Working Group (HI-TWG) for their oversight role during the preparation of this manual.

A handwritten signature in black ink, appearing to read 'Henry Mwebesa', written over a horizontal line that tapers to a point on the right. There are some small dots below the line.

Dr. Henry Mwebesa
Director General for Health Services

Acronyms

5S	Five steps of Sort, Set, Shine, Standardize and Sustain
AC	Alternating Current
AEO	Assistant Engineering Officer
BP	Blood Pressure
BME	Biomedical Engineer
CD4	Cluster of Differentiation 4
CHS(HID)	Commissioner for Health Services, (Health Infrastructure Department)
CO ₂	Carbon Dioxide
CT	Computerized Tomography
CW	Central Medical Equipment Maintenance Workshop
DC	Direct Current
DHOs	District Health Officers
ECG	Electrocardiogram
ENT	Ear, Nose and Throat
HA	Hospital Administrator
PHA	Principal Hospital Administrator
HC	Health Centre
HGB	Haemoglobin
HID	Health Infrastructure Department
HI_TWG	Health Infrastructure Technical Working Group
HIV	Human Immunodeficiency Virus
HMIS	Health Management Information System
JICA	Japan International Cooperation Agency
LCD	Liquid Crystal Display
LED	Light Emitting Diode
ME	Medical Equipment
MoFPED	Ministry of Finance, Planning and Economic Development
MOH	Ministry of Health

MRI	Magnetic Resonance Imaging
N ₂ O	Nitrous Oxide
NACME	National Advisory Committee on Medical Equipment
O ₂	Oxygen
OPD	Out Patient Department
OT	Operation Theatre
PDU	Procurement and Disposal Unit
PMT	Photo Multiplier Tube
PNFP	Private–Not–For–Profit
PPDA	Public Procurement and Disposal of Public Assets Authority
PCM	Planned Corrective Maintenance
PPM	Planned Preventive Maintenance
QI	Quality Improvement
RBC	Red Blood Cell
RF	Radio Frequency
RPM	Rotations Per Minute
RRHs	Regional Referral Hospitals
RWMC	Regional Workshop Management Committee
RW	Regional Medical Equipment Maintenance Workshop
SOPs	Standard Operating Procedures
SpO ₂	Saturation of Peripheral Oxygen (Pulse Oximeter Oxygen Saturation)
SVR	Sample Rotary Valve
WBC	White Blood Cell
WIT	Work Improvement Team

Chapter 1: Introduction

1.1 Background

Medical equipment plays an important role in our Health Care system and there are numerous equipment adapted to Health Care. In the health facilities, medical equipment will range from injection needles, centrifuges, stethoscopes, blood pressure machines, to X-Ray Machines, MRI etc. designed to help medical personnel carry out diagnosis and treatment of patients.

Optimal performance of medical equipment is required to ensure safety, accuracy and expected results. To keep medical equipment in good working condition and optimal performance is the function of a medical equipment maintenance unit in a health facility.

The Ministry of Health (MOH) set up Regional Medical Equipment Maintenance Workshops (RWs) to carry out medical equipment maintenance in public health facilities. For Private–Not–For–Profit (PNFP) facilities, the workshops will provide maintenance support, whenever necessary and they shall be required to procure the necessary spare parts. In the hospitals, there are different engineering personnel deployed to man hospital maintenance units and carry out maintenance on medical equipment.

This operation manual is intended to support and guide engineering personnel in the RWs and Hospital maintenance units in the maintenance of essential medical equipment. It will also guide the equipment users and technicians on the operation and working principles of some commonly used equipment.

The second edition of the operation manual for RWs has been prepared in two volumes, namely:

- Volume I: Operation Manual for Regional Medical Equipment Maintenance Workshops and Medical Equipment Maintenance Guidelines
- Volume II: Standard operating procedures for carrying out PPM on selected commonly used equipment.

1.2 Overall Objective

The overall objective of preparing this manual is to define the roles of RWs and other stakeholders; and provide guidelines for medical equipment maintenance management including maintenance planning, budgeting, equipment user training and disposal of medical equipment.

1.3 Specific Objectives

The specific objectives of this manual are to:

- Define the role and function of RWs and Central Medical Equipment Maintenance Workshop (CW)
- Identify and define the roles of various stakeholders
- Streamline the operations of RWs
- Provide maintenance guidelines for commonly used medical equipment in health care facilities
- Assure safety, proper use, care and management of commonly used medical equipment
- Provide Standard Operating Procedures (SOPs) for carrying out Planned Preventive Maintenance (PPM) on commonly used medical equipment in health care facilities **attached as Volume II to this manual.**

1.4 Stakeholders in Medical Equipment Management

The following categories of health care workers and providers play a significant role in the management of medical equipment.

- **Users of Medical Equipment**

Equipment Users includes, doctors, nurses and other allied health professionals whose primary function is to use medical equipment for diagnosis and treatment of patients. Their main role will be to care for the equipment including reporting equipment failure to the maintenance unit.

- **Medical Equipment Maintenance Personnel**

Medical equipment maintenance personnel include, artisans, technicians and engineers employed by health facilities or MOH to offer equipment maintenance services. Their main role is to identify and isolate the fault, and take remedial action to restore full functionality of the medical equipment.

- **Hospital Administrative Personnel**

This category includes the hospital administrators, procurement personnel, accounts personnel and other staff responsible for safe custody of hospital stores. This category is involved in planning, assets/stores management, procurement of spares, and financial management.

- **Medical equipment manufacturers and vendors**

The role of this category is to offer spare parts for the equipment they manufacture or sell and offer After Sales

Services including contract maintenance services for specialized equipment.

- **Development and Implementing partners**

They support in the area of funding, carry out capacity building to maintain medical equipment

- **Local governments**

They fund and supervise medical equipment maintenance, management and disposal.

1.5 Definition of Medical Equipment Maintenance

Medical equipment maintenance refers to a set of activities conducted on a medical equipment to keep it in optimum working. It is divided into two major categories namely:

- **Planned Preventive Maintenance (PPM)**

This refers to regular safety and performance inspection carried out on equipment to evaluate risk and reduce failure so as to enhance its performance, efficiency and reliability.

It involves cleaning, regular function/safety tests and making sure that any problems are picked up before they cause a breakdown.

PPM is recommended for all medical equipment. It will enhance the efficiency, effectiveness and reliability of medical equipment and must be carried out at appropriate frequency as suggested by the manufacturer or based on workload and usage.

- **Corrective Maintenance (or repair)**

Corrective maintenance is a task performed to identify and rectify a fault on a broken-down equipment, machine or system to restore to it to its original operational condition.

The choice of approach for preventive and corrective maintenance depends on the complexity of equipment, equipment usage and/or cost benefit analysis by the in-house trained engineering personnel. The majority of equipment in our health care system is basic in nature and can be handled in-house by any engineer/technician at RW level with training in medical equipment maintenance.

For specialized and advanced equipment, the vendor should provide maintenance services through a combination of on-call services and a maintenance contract negotiated at the time of purchase. Maintenance contracting should be reserved for medical equipment for which there is no in-house capacity to maintain or when specialized tools and technical expertise are needed and it's not economically viable to develop this capacity in-house.

- **Planned Corrective Maintenance (PCM)**

This refers to corrective maintenance that is planned and performed to rectify a fault on a broken-down equipment or system to restore to it to its original operational condition.

Chapter 2: Organization Structure for Medical Equipment Maintenance

There is a wide range of medical equipment and hospital plants at different levels of the health care delivery system. The staffs in these health facilities are responsible for ensuring that it is used and stored properly; and faulty equipment is reported to departments responsible for maintenance.

At the central level, the Health Infrastructure Department (HID) under the Directorate of Policy, Strategy and Development is responsible supervising and monitoring equipment maintenance in health facilities; and providing policy guidelines on operation and management of medical equipment. On the other hand, the National Advisory Committee on Medical Equipment (NACME) is mandated to give appropriate policy advice to MoH regarding medical equipment.

At District and Health facility levels, the management of medical equipment is the responsibility of the respective Administrative and Technical Departments. The District Health Officers (DHOs) are responsible for planning and budgeting for management of medical equipment in the District Local Governments.

The RWs under the Regional Referral Hospitals (RRHs) are mandated to carry out medical equipment maintenance in health facilities in their catchment areas. However, the NRHs and Specialized Hospitals have fully fledged Engineering Departments with Engineers and Technicians who maintain their equipment.

On a case by case basis, the RWs, NRHs and Specialized Hospitals may seek support from the Central Medical Equipment Maintenance Workshop (CW) as a referral workshop.

2.1 National Advisory Committee on Medical Equipment (NACME)

The main function of NACME is to guide the MoH on the country's medical equipment needs and develop appropriate policy framework.

This includes advising on equipment specifications, acquisition methods, equipment standardization and maintenance of medical equipment; bearing in mind efficiency, cost-effectiveness and appropriateness of technology.

2.2 Health Infrastructure Department (HID)

The HID has two major Divisions with the mandate to formulate policies, strategies and guidelines on health infrastructure development and management.

The Divisions are:

- Civil and Sanitary Engineering
- Biomedical and Electro-Mechanical Engineering

1) Civil and Sanitary Engineering Division

All Civil and Sanitary Engineering services in the Sector are coordinated by this Division. It provides technical support and supervision of equipment pre-installation works and ensures that fixed equipment is installed safely.

2) Biomedical and Electro-Mechanical Engineering Division

This Division encompasses Electrical, Mechanical and Biomedical engineering disciplines. The Electrical, Mechanical and Biomedical engineering professionals are responsible for preparing appropriate equipment

specifications and ensuring that procured equipment conforms to national and international standards; and that the equipment is maintainable.

The key equipment related activities handled by this Division include:

- i) Supervision and monitoring the management of the complete life cycle of medical equipment, hospital plants and medical furniture in public health facilities.
- ii) Liaising with NACME to update the national medical equipment policy, standard equipment lists and technical specifications.
- iii) Coordinating equipment inventory taking, update and quantification of needs of health facilities.
- iv) Organizing capacity building training for engineers and technicians to keep abreast with the fast-changing trends in healthcare technologies.
- v) Co-ordination and management of service contracts.
- vi) Planning and budgeting for RWs capacity building needs
- vii) Assessing the performance of RWs.
- viii) Coordinating development and implementing partners.

2.3 Regional Medical Equipment Maintenance Workshops

In order to decentralise medical equipment maintenance, the MoH established Regional Medical Equipment Maintenance Workshops (RWs) at RRHs to maintain medical equipment in health facilities under their catchment area.

To date there are twelve (12) RWs located at Arua, Jinja, Lira, Gulu, Soroti, Mbale, Hoima, Fort Portal, Kabale, Mubende, Moroto, Masaka RRHs and Central Medical Equipment Maintenance Workshop in Kampala. More RWs are to be established in RRHs that do not have regional workshop.

The Central Medical Equipment Maintenance Workshop operates under the stewardship of HID and is a referral workshop for all other RWs. It supports other RWs, coordinates national training activities and supervises centrally managed maintenance contracts-




Fig 1: Medical Equipment Maintenance Organizational Framework




Fig 2: Location and Catchment Areas of RWs

2.4 National Referral Hospitals (NRHs)

At the National Referral Hospital level, there is a fully-fledged Engineering Department that works independent of the RWs. It is responsible for health infrastructure planning, budgeting and management for NRHs.

2.5 Regional Referral Hospitals (RRHs)

At RRHs, the management of medical equipment and hospital plants is the responsibility of the Hospital Administrators and respective RW Manager. While the RWs are part of the Hospitals, their main mandate is to maintain equipment in the RRH and health facilities in their catchment area. This is one of the outreach services of the RRHs.

2.6 District Health Services

At the District level, the management of medical equipment is the responsibility of the DHO and the respective in-charges of the Health Facilities. District Engineers in Local Governments assist the DHO to plan for health infrastructure development and management.

Biomedical Engineering services are provided by the team from the RWs to carry out periodical equipment maintenance in the General Hospitals, HCIVs, HCIIIs and HCIIIs. The facility-based Engineers/Technicians are responsible for carrying out basic health infrastructure maintenance using the facility budget.

Chapter 3: Function, Operation and Role of RWs

The CW and RWs were established to carry out medical equipment maintenance in health facilities as mandated in the national medical equipment policy.

The RWs are a support service section under the RRH to maintain equipment in the lower level health facilities as an outreach service of the RRH. A RW set up consists of;

- A workshop building with office space, stores and work area – **See Annex 1 for the standard floor plan of a RW**
- A purpose-built mobile workshop vehicle or van
- Tools, test equipment and office furniture – **See Annex – 2 for list of tools, test equipment and furniture**
- Staff consisting of Engineers, Assistant Engineering Officers (AEO), Engineering Technicians and Artisans.




Fig 3: RWs Macro Organogram

Operationally, Engineers and Technicians from the RW visit Hospitals, DHO's stores and HCIVs to carry out equipment maintenance on site. For HCII and HCIII, the faulty medical equipment is to be delivered to the nearest hospital, DHO's stores or HCIV by the In-charge of the HC for the workshop team to carry out maintenance. On a case by case basis however, the RW team may visit a HCII or HCIII to repair immovable equipment like a generator, solar system. Additionally, all HCIII laboratory hubs and high volume HCIII with electrical equipment should be visited by the RWs maintenance teams to maintain the equipment on site.

3.1 Central Medical Equipment Maintenance Workshop (CW)

The CW in Kampala is the RW for health facilities in the central region and a referral workshop for work that cannot be handled by Engineers/Technicians from other RWs. It is also a training centre for hospital-based technicians. Other roles and functions of the CW include the following:

- i) Supervision of all RWs through the Biomedical and Electro-Mechanical Engineering Division of HID focusing on maintenance planning and execution of maintenance activities – i.e. review of work plans and budgets, maintenance schedules and quarterly reports.
- ii) Identification of suitable service providers and guidance on appropriate procurement methods for medical equipment spare parts.
- iii) Plan and coordinate training of RW staff and mentorship.
- iv) Support other RWs to carry out and update medical equipment inventory.
- v) Provide specialized medical equipment testing and calibration services. (e.g. mass, temperature, rotation speed, pressure and electrical safety)

3.2 Regional Medical Equipment Maintenance Workshops

The RWs are based at the RRH and their main mandate is to carry out maintenance of medical equipment in health facilities within the catchment area of the RRH.

Each RW shall be headed by a Workshop Manager designated by the Hospital Director in consultation with the CHS(HID). The Workshop Manager shall be an Engineer. He/she will be deputized by an Engineer or AEO with good knowledge of medical equipment maintenance and computer skills. He/she will be responsible for medical equipment inventory management and user training.

While day to day supervision of the operations of RWs is under the Administration of the RRH, each workshop shall have a RW Management Committee (RWMC) that over sees its operations.

The RWMC will consist of the HD/RRH, RW Manager, Chief Administrative Officers (CAOs), DHOs, Hospital Medical Superintendents (MSs), Heads of Health Sub-Districts and the District Local government Secretaries for Health.

The main functions of RWs include the following:

- i) Maintenance of medical equipment in health facilities in RRH/RW catchment area.
- ii) Supervise and/or carry out installation of new equipment in health facilities in the RRH/RW catchment area.
- iii) Collect and update medical equipment inventory in health facilities in their region.
- iv) Advise health facility managers on medical equipment status, equipping gap and disposal of obsolete equipment.

- v) Carry out equipment user training (proper use, operation, handling and first line maintenance).
- vi) Preparation and submission of quarterly RW performance reports to the HID/MoH.
- vi) Organizing and participating in the RWMC Meetings to share with stakeholders the RW performance and challenges.

3.3 RW Maintenance Work Flow Chart

To ensure effective monitoring of workshop activities and establish an efficient maintenance system, the work flow described in the chart below will be followed. For this to be effective, the Users/Health facilities must play their part by promptly reporting equipment breakdown.




Fig 4: RW Maintenance Work Flow Chart

Chapter 4: Planning, Budgeting and Management of Medical Equipment Maintenance

Lack of medical equipment maintenance results in equipment breakdown and prolonged downtime affects the reliability and quality of health care. Lack of a maintenance policy can result in poor planning, lack of maintenance budgets and thus no spares parts and accessories. Many health facilities have faulty equipment because the operation and maintenance (O&M) requirements are not planned for in advance. This renders many equipment unusable and several devices lie idle because of lack of consumables, spares parts or inadequate funds for O&M.

The staff of the RWs are responsible for planning and budgeting for maintenance activities. For day to day duties, the RW staff are responsible to the Hospital Director and Hospital Administrator. At Ministry level they are supervised by Engineers in the HID/CW, in Kampala.

The following activities will promote effective maintenance planning and management:

- i) Regular medical equipment inventory taking, update and analysis of data.
- ii) Use of updated and analysed data for maintenance planning, budgeting and preparation of maintenance schedules.
- iii) Preparation of quarterly and annual work plans and budgets.
- iv) Preparation of maintenance schedules and implementation of PPM.
- v) Preparation of maintenance reports including job cards, service stickers, quarterly and annual performance reports.

- vi) Holding regular meetings – i.e. workshop staff meetings, regional workshop steering committee meetings and RWs performance review meetings.

4.1 Medical Equipment Inventory Collection and Update

The Workshop manager shall endeavor to collect and update inventory of all medical equipment in each public health facility in the catchment area of the RW. Proper entry should be made in the New Order for Managing Anything Data (NOMAD) inventory database or any other software approved by the HID/MoH. The inventory record should contain the details in the Annex-4:

- i) Name of health facility, department, section and room
- ii) Equipment Name and Model Number
- iii) Manufacturer (with contact details where possible)
- iv) Serial No. of equipment or allocated identification No.
- vi) Date/Year of Installation/purchase/manufacture
- vii) Equipment functional status “Condition A–F” – **see Annex 3**.

Note: During inventory taking, medical equipment shall be categorized into the following functional conditions:

- A:** Good and in use
- B:** Good but not in use
- C:** In use but needs repair
- D:** In use but needs replacement
- E:** Out of order but repairable
- F:** Out of order and should be replaced

In order to have up to date inventory data for maintenance planning, medical equipment inventory taking and update should

be carried out on a continuous basis in accordance with the mandatory annual inventory update.

Equipment inventory update should be done whenever new equipment has been delivered in a health facility, the equipment working condition has changed or equipment has been transferred from one location to another.

4.1.1 Inventory taking Rules and Procedure for Inventory Update

An accurate and up to date inventory database is important for proper maintenance planning and scheduling. It should be carried out in good time to inform decision making on equipping gaps, disposal, maintenance planning and budgeting. The following methodology shall be followed when carrying out inventory:

- Visit each department and section of the hospital (**including the stores**) and record all available equipment against the NOMAD generated inventory list for the health facility and/or the user's inventory record if available.
- Carry out a physical count and inquire whether the quantity is adequate for the current staffing and workload.
- Consult the user on the current functional status (i.e. whether the equipment is working or not) and use.
- Collect/verify and record all required information as detailed in the Medical Equipment Inventory Form – **See Annex 4**

The RW inventory data collection team in consultation with the hospital staff should identify equipment that is due for decommissioning and/or disposal.

4.1.2 Basic Medical Equipment Inventory Data Analysis and Use

Proper maintenance planning by RWs can only be achieved using analysed information from collected equipment inventory data. In

order to efficiently utilize the limited funding availed to RWs, evidence-based planning, budgeting and maintenance scheduling is recommended using up-to-date inventory data.

All RWs shall carry out the following inventory data analysis to enable use of the data for efficient and cost-effective maintenance planning and budgeting.

Data analysis based on equipment functional condition	Application of information for planning
Category “A”	<ul style="list-style-type: none"> • Used to identify all equipment that needs PPM and their location • Used to prepare a PPM schedule and budget – i.e. time required to carry out PPM, required spare parts & consumables
Category “B”	<ul style="list-style-type: none"> • Used to identify all idle equipment in good working condition and the reasons for not being in use. • Used to identify all idle equipment for which user training is required. • Used to prepare a user training plan and execution budget.
Category “C” and “E”	<ul style="list-style-type: none"> • Used to identify all equipment that needs repair and their location. • Used to prepare the list of spare parts and supplies required to fix faulty equipment. • Used to prepare a PCM schedule and budget.
Category “D”	<ul style="list-style-type: none"> • Used to identify old, inefficient and expensive equipment to maintain; and their location.

Data analysis based on equipment functional condition	Application of information for planning
	<ul style="list-style-type: none"> • Used to prepare a list of equipment that should be replaced in the short term.
Category “F”	<ul style="list-style-type: none"> • Used to prepare lists of equipment that should be disposed of and replaced.

4.2 Work Plan and Budgeting for Medical Equipment Maintenance

In order to plan for maintenance activities, the Workshop Manager should prepare a work plan and budget. The work plan should have an itemized summary of planned activities, targets, timelines and intended expenditure on a quarterly basis. The work plan and budget should indicate the proposed source of funding.

Adequate time should be allocated for procurement of goods and services including spare parts. The standard templates for preparation of work plans should be used – **see Annex 5 for Work Plan Preparation template.**

Currently, medical equipment is maintained using resources from the central treasury through RRH and Local Governments. Additional funding for equipment maintenance is provided for by Development & Implementing partners.

For regional medical equipment maintenance in the catchment area for the RRH, workshop funds are budgeted under Program 03 (Regional Maintenance). This is a ringfenced budget for medical equipment maintenance for the whole region and not just the RRH. It should not be used for civil works in the RRH except for minor renovation works on the RW. This budget excludes the wages of the RW staff.

4.2.1 Process of Preparation of Annual Work Plans and Budgets

The budgeting process starts in October of every year when the MoFPED sends out a budget call circular.

Each RRH makes an annual work plan and the RW workshop plan and budget is part of it. The RRH submits the work plan and budget to MoH for onward submission to MoFPED which makes budgetary allocation for various Ministries and/or Departments.

During this time, the Workshop Manager should prepare and submit the RW priority activities, proposed budget allocations and implementation schedule; to the Hospital Director for incorporation in the overall work plan for the RRH.

4.2.2 Preparation of RW Annual and Quarterly Budget

Since funds are disbursed to RWs on a quarterly basis; planned activities should be scheduled across the year in a logical manner. Identify routine activities that are repetitive and spread them across each quarter and add other activities that can be carried out in the available time and allocated funds for each quarter. While it is usually normal to divide the annual budget allocation by 4, it is important to critically look at the processes and inputs (i.e. time, human resources and finances) needed to realize planned activities. The RW budget should provide for the following:

- i) Procurement of spare parts and maintenance workshop supplies
- ii) Labour cost including outsourced services
- iii) Mobile workshop operational costs (i.e. fuel, tyres, servicing and repair)

- iv) Subsistence costs while on maintenance trips (i.e. Per diem , Safari day allowances)
- v) RWMC Meetings
- vi) RWs performance review Meetings
- vii) Training for users and workshop technicians
- viii) Replacement of essential tools and consumables
- ix) Office stationery, ICT services and supply.
- x) Maintenance of office equipment
- xi) Payment for telephone, water and electricity bills.

4.2.3 Preparation of Quarterly Maintenance Schedules

As one of the planning tools, each RW will prepare a quarterly maintenance schedule. The maintenance schedule will identify the Health facilities to be visited, allocated time and human resources.

The quarterly schedule shall be prepared after confirming the available funds and prioritized planned activities for the quarter.

The RW manager shall send out the maintenance schedule in good time to allow the Health facility managers plan and mobilize the staff to avail the equipment that needs maintenance.

The maintenance schedule should be followed as much as possible by the RW to ensure efficiency and cost effectiveness.

4.2.4 Budgeting for Emergency Works

There will always be cases of critical equipment breaking down and the RW may be expected to respond immediately. To cater for such incidences, a contingency (10% to15% of the overall

quarterly budget allocation) should be provided for in the quarterly work plan to cater for emergency callouts. The contingency funds should be spent on a case by case basis on express approval of the RRH Director; and could be drawn from the petty cash and replenished immediately.

4.3 Medical Equipment Maintenance Planning and Management

In the management of RW operations and Equipment Maintenance, it is important to involve stakeholders especially the Users and Administrators. This will help each player to have a good understanding of WHAT to do, WHEN to do it, WHO is to do what; and to agree on the necessary MEANS to do it.

The following scheduled reports and meetings will provide opportunities for engagement and sharing experiences and challenges of medical equipment management:

- Quarterly performance reports
- RWs performance review meetings
- RWMC meetings
- RW Staff meetings

RWs will implement 5S activities designed to improve the work environment, safety, and work flow. Additionally, RWs will implement Continuous Quality Improvement (CQI)–Total Quality Management (TQM) activities aimed at continuously improving work processes and solving identified bottlenecks within the means of the workshop human resources and budget. This Manual has a summarized description of 5S implementation steps for the workshop environment and an Action Matrix for implementing small CQI activities attached at **Annex 6 (a) and Annex 6(b)**. The detailed description of 5S–CQI–TQM implementation framework can be found in the 5S–CQI–TQM Hand Book and Implementation Guidelines produced by MoH.

The key planning tools for RWs shall be;

- 1) Annual and quarterly work plan and budget
- 2) PPM schedules
- 3) CQI Action matrix

4.3.1 Preparation of Quarterly Reports

The Workshop Manager shall prepare and submit quarterly performance reports to the CHS(HID) through the Hospital Director every 15th day of the months of April, July, October and January of each year.

The quarterly report shall highlight the overall progress in implementation of planned activities, set targets and performance stating the achievements, outputs and funds used to obtain the outputs. The quarterly report will contain the following:

- Executive summary – highlighting outputs/outcomes, analysis of equipment working condition, budget performance, challenges and recommendations to improve performance.
- Planned activities, set output targets and status of implementation.
- Summary of medical equipment maintained in each health facility.
- Expenditure during the quarter clearly separating expenditure on spare parts, per diem, workshop operational costs (i.e. electricity, water, communication as necessary), labour, and vehicle maintenance,
- Planned activities and expenditure for the next quarter.
- Challenges met
- Recommendations
- Conclusion

A standard format for preparing quarterly reports is attached as Annex 7 and attachment annexes 7(a), 7(b), 7(c), 7(d) and 7(e).

All Job Cards for work carried out in the standard format should be attached – See Annex 8

4.3.2 RWs Performance Review Meetings

The Managers of RWs and their RRH Supervisors will meet at the end of every quarter to discuss issues relating to Medical Equipment Maintenance and Management of RWs.

- **Purpose of meeting**

The meeting is for reviewing performance of the RWs, sharing experiences, challenges and charting a way forward to improve maintenance of medical equipment across the country.

- **Frequency**

The CHS(HID) through the CW shall coordinate RWs to organize and hold quarterly performance review meetings hosted by a RW agreed upon by all the RW on alternate basis.

- **Participants**

The RWs performance review meeting will include Workshop Managers and RW staff, HD/PHA of RRHs, Senior staff of HID/CW and Development Partners supporting equipment maintenance; and other stakeholders as shall be deemed fit from time to time (e.g. user trainers, representatives of manufacturers, equipment vendors, Service providers).

Each RW manager shall prepare a Power Point presentation and share the summary of the quarterly performance. Minutes of the meeting shall be prepared by the host RW.

4.3.3 Regional Workshop Management Committee (RWMC)

In order to bring together stakeholders involved in equipment use and management, each RW shall hold a RWMC meeting at least twice a year.

- **Purpose**

The RWMC meeting will be to review the performance of the RW in respect to equipment maintenance, share experiences, challenges and chart a way forward to improve performance of the RW. The work plan, budget and maintenance schedules will be reviewed and approved for each year and/or quarter. The meeting will be chaired by host District/Hospital while the RW will be the secretariat for the meeting.

- **Frequency**

The RWMC meeting shall take place **at least twice a year** to coincide with the review of the quarterly/annual report prepared by the RW.

- **Recommended participants**

The following stakeholders will constitute the RWMC membership:

- HD/PHA of the RRH
- Medical Superintendents/Hospital Administrators of the beneficiary Hospitals
- DHOs of the beneficiary Districts
- 25% of CAOs and Heads of Health Sub-Districts on a rotational basis
- District Local government authorities responsible for health
- Workshop Manager and Regional User Trainer
- Representative of HID/MOH
- Development Partners supporting equipment maintenance

Note: The meeting should always review the proposed work plan, schedule and budget for next quarter/year.

4.3.4 Internal Workshop Staff Meetings at workplace

RWs are manned by a Workshop Manager and other technicians in the RRH. To ensure cohesion and coordination in the RW, the workshop manager will hold regular internal meetings with other staff. The RW Manager will use these meetings to discuss action plans, allocate work and review progress and challenges.

- **Frequency**
Internal meetings between the workshop manager and his /her staff should held on a weekly basis to ensure timely follow up of assigned tasks to be performed on a daily, weekly and/or monthly basis.
- **Discussion points**
Relevant tasks, teamwork, staff welfare, assessing progress of CQI activities, challenges, performance improvement and agreeing on timelines by all staff members.
- **Meeting Minutes and Work Reports**
For every meeting, there shall be minutes of the meeting. The minutes will document Attendance, issues discussed and agreed positions, and a summary of follow up Actions and responsible persons.

The Workshop manager will put in place a mechanism to ensure that action points agreed upon in the internal meetings are documented and followed up.

For every maintenance work carried out, the Engineer/Technician shall prepare a report documenting the tasks he/she undertook for the day/week and the overall status of the task.

Minutes of all meetings held shall be put on file for future reference. **Chairing of meetings and making of Minutes shall be a shared responsibility by all staff of the workshop on a**

rotational basis including the workshop manager as a strategy to develop the capacity of all staff.

Note:

1. Daily meetings should last not more than 30 Minutes.
2. Meetings held on Weekly basis should not last more than 1 hour.

Chapter 5: Procurement and Disposal of Goods and Services by RWs

While RWs will get involved in procurement and disposal of medical equipment; this is not their core function. This is the function of the respective Heads of the Health facilities and the Procurement and Disposal Units (PDU). Their role in medical equipment procurement should be limited to;

- Review and/or provision of technical specifications
- Inspection and testing equipment to confirm conformity to contract technical specifications and functionality
- Installation and commissioning new equipment
- Recommending medical equipment for decommissioning and disposal based on technical considerations (e.g. age, obsolescence, lack of spare parts, reliability and safety).

For day to day workshop operations, the RWs will be involved in procurement of medical equipment spare parts and maintenance services. In this case, the workshop staff will be involved in procurement processes as a User Department.

5.1 Management of Procurement by RWs

Procurement of spare parts and services by RWs will be guided by the Public Procurement and Disposal of Public Assets Authority (PPDA), the Regulation thereto and guidelines in force. All RW managers need to acquaint themselves with PPDA guidelines and the different procurement methods. RWs will be responsible for planning for spare parts needs (i.e. quantification, specifications and budget), out sourcing maintenance services on a case by case basis and disposal of spare parts that are not useful.

As a User Department, the RW will play the following procurement roles;

- i) Determine and quantify workshop requirements
- ii) Propose/provide technical specifications for the required supplies
- iii) Raise procurement requisitions using the Procurement Form 5 (Requisition for approval of procurement).
- iv) Participate in the selection of supplies and service providers
- v) Participate in inspection, testing and verification of purchased supplies.

Note:

1. RW staff must exercise high moral and ethical conduct while handling procurement of goods and services
2. No supplies should be used before they are taken on charge in the stores and issued out.

5.1.1 Determination of Spare Parts Needs

Determination of spare parts needs shall be based on demonstrated maintenance needs and requirements for Planned Preventive Maintenance (PPM). Annex 3 of the RW quarterly performance report will always be useful in determining spare parts for pending Jobs and drawing up a PCM schedule.

Problem oriented and evidenced based planning and budgeting will be the basis for maintenance planning and scheduling to ensure efficient and cost-effective utilization of available resources.

- Medical equipment inventory for each health facility will be analyzed to determine the maintenance condition of the equipment. Repairs carried out will always be supported with a **Job Card – see Annex 8**.

- RWs will endeavor to print serialized Job Cards and enforce the preparation of Equipment breakdown Report Forms (**HMIS Form 011**) by health facilities.
- A list of fast moving spare parts will be prepared and used to determine quarterly requirements. Where possible a list of fast-moving spare parts for each equipment will be developed from time to time.

5.1.2 Procedure for Purchasing Spare Parts

Procurement of spare parts and tools shall be carried out in line with the PPDA guidelines. The RW will identify and quantify the spare parts requirements, and raise the Procurement requisition to start the procurement process through the PDU.

While the RW will allow the PDU to manage the procurement process, it is advisable for the RW manager to assist the PDU to identify possible suppliers and appropriate procurement method. This is important because of the specialized nature of medical equipment spare parts.

- Some spare parts are manufacturer specific and direct procurement would be the most appropriate procurement method.
- Also, prequalification of suppliers would help reduce the time needed to identify and place orders.
- Framework contracts should be signed for supply of fast moving medical equipment spare parts to minimize time spent in procurement processes and facilitate fast placement of orders. This should be the preferred procurement method for spare parts purchase.

Whatever procurement method is used, it is important that the process is transparent, fair, efficient and cost effective.

5.1.3 Stores Management by RWs

RWs have stores and must maintain inventory records. Management of all stores (i.e. tools and spare parts) shall be in accordance with standard stores management practices. The RW shall secure and maintain stores “Stock cards”, “Stores Requisition Forms” and “Stores Issue Forms”.

Stock cards shall be maintained weekly, while monthly stock taking will be carried out for all inventories. No spare part shall be issued out without clearly identifying it (i.e. name, part number and model).

5.2 Disposal of Items

Disposal is part of the procurement process and must be planned for. The life cycle of any equipment ends with disposal; but it is often difficult for users to decide when to decommission and how to dispose of equipment.

The reasons for decommissioning and disposal of equipment will include:

- Wear and tear beyond economic repair
- Damage beyond economic repair
- Unreliability
- Safety concerns
- Clinical or technical obsolescence
- Unavailability of spare parts
- Availability of more cost-effective equipment or clinical procedure or technology

Once items for disposal are identified, the PPDA guidelines and procedure should be followed. The disposal process would include;

- A) Identification of items to be disposed of
- B) Submission of the list of items for approval by the Accounting Officer (i.e. Hospital Director and/or Hospital Board).
- C) Establishment of a Board of Survey team to verify and approve the list of items for disposal.
- D) Initiation of disposal process through preparation and submission of Form 28 (Requisition for Disposal) to the PDU to start the disposal process.

Note: The role of the RW will be advisory as the Accounting Officer/PDU has the mandate to dispose of assets/ stores.

Chapter 6: Support Supervision, Monitoring and Evaluation of RWs

The Quality Assurance Department of MoH is responsible for monitoring and evaluation of healthcare services. This is done through the Area Teams which assess and monitor implementation of various policy guidelines, planned activities and performance against Sector indicators.

The Evaluation and Monitoring assessment form for Area Teams includes Health Infrastructure and RW staffs need to provide the information on medical equipment maintenance and its impact on healthcare delivery.

6.1 Support Supervision for RWs

For the RWs, specific technical support supervision, monitoring and evaluation will be conducted by the HID/MoH and focus will be put on the following;

- Maintenance planning and budgeting. Evidence based planning and budgeting should be the basis for budget allocation.
- Availability of adequate spare parts. Spare parts stock planning and control should ensure that essential medical equipment is well maintained at all times.
- Efficiency and cost effectiveness. Operational costs should be optimized to a minimum level but with high outputs and outcomes.
- Prudent financial management and planning. Maintenance scheduling and resource allocations should be optimized to ensure efficiency and cost effectiveness.
- Equipment downtime and response time. The response time to emergency calls for repair of life support equipment should not be more than 24 hours. The PPM schedule

should aim at keeping at least 70% of the available medical equipment in good maintenance condition.

1) Frequency of visits

The HID will carry out support supervision and monitoring of RWs every quarter.

A standard supervision, monitoring and RWs assessment sheet has been designed to assess RWs – **see Annex 9**.

Chapter 7: Guidelines for Medical Equipment Maintenance

7.1 Introduction

The main objective for any maintenance system is to ensure prolonged use of available equipment over its design life to provide safe and reliable healthcare.

All RWs shall ensure proper maintenance of medical equipment to sustain the intended healthcare benefits and to preserve capital investments. To achieve this, medical equipment shall be maintained in working order and periodically calibrated to ensure accurate results.

To ensure efficient maintenance of equipment, RWs will plan and budget for maintenance under two main categories:

- Planned Preventive Maintenance
- Corrective maintenance

7.2 Planned Preventive Maintenance (PPM)

PPM is usually scheduled at specific intervals and includes specific maintenance activities such as lubrication, calibration, cleaning of filters or replacement of spare parts that are expected to wear out after stipulated time or workload (e.g. bearings, tubings).

All RWs shall endeavor to implement PPM and develop maintenance schedules based on the principle of “Problem oriented maintenance planning and budgeting” to optimize utilization of the limited funds and human resources.

Effective planning for preventive maintenance will involve careful selection of the equipment to be included in the plan. Evidence based decisions must be made while deciding equipment to include in the maintenance schedule in order to reduce costs. Maintenance planning shall always ensure that essential medical equipment for basic diagnosis, infection control, surgery and treatment are kept in good working condition. The overriding considerations will always be safety, efficiency and cost effectiveness.

RWs will carry out equipment calibration and testing using the SOPs developed by the maintenance team. Proper documentation of the testing/calibration results shall be recorded and filed together with the corresponding Job card issued out for the job.

7.3 Corrective Maintenance Services

Corrective maintenance refers to corrective actions undertaken in the event of breakdown of a piece of equipment. In this case, the equipment is repaired or calibrated after failure.

At all times, medical equipment in use should be free from any fault regardless of how minor the fault is and all repair work should be performed by a competent engineer/technician.

The user department should:

- 1) Record details of the defect
- 2) Fill-in a “Maintenance Requisition Form”. Refer to **Annex 10**.
- 3) Contact the RW engineer/technician who should in turn decide whether to carry out the repairs or to contact the maintenance contractor.

- 4) Ensure that information regarding equipment breakdown is passed to all staff during shift change and the head of the institution.

All equipment breakdown occurring on the wards should be recorded on the Maintenance Requisition Form mentioned above.

7.4 Maintenance Policy Guideline

The choice of implementing Preventive or Corrective Maintenance depends on the complexity of equipment, availability of After Sales service and cost. However, it is often said that “prevention is better than cure”.

Maintenance of medical equipment by the RWs will be based on PPM for all essential medical equipment that can be handled by RWs; and a mix of PPM and Corrective maintenance for medical equipment that need outsourcing the maintenance services.

Equipment maintenance shall be carried out by both In-House trained engineers/technicians and the Manufacturer’s representatives or Vendor.

For specialized and advanced equipment, the HID/MoH and RWs will on a case by case basis outsource maintenance services through a combination of on-call services and maintenance contracts negotiated and signed with competent service providers.

Where maintenance contracts exist; the RW Manager will be required to supervise, monitor and to the extent that may be practical, certify work done by the Engineers/Technicians of the service providers especially for maintenance in the RRHs.

7.5 Levels of Maintenance

For purposes of this Manual, three levels of maintenance will be observed:

Level 1: First-line by Equipment User

This refers to maintenance activities that can be carried out by the user. This will include dusting equipment, cleaning filters, and checking power supply source without opening the unit and without moving it away from the point of use.

Other first line maintenance activities could include the following:

- Equipment decontamination and sterilization
- Functional checks
- Calibration checks
- Safety checks

Note: First line maintenance may be carried out by Hospital based Technicians too.

Level 2: By Engineers and Technicians

This refers to maintenance carried out by an engineer/ technician when first-line maintenance cannot rectify a fault or when a regular scheduled check, PPM and/or calibration are due.

Level 3: Specialized Maintenance by Technician/Engineer or Manufacturer's representative/Vendor

This refers to maintenance activities that need higher level technical expertise, troubleshooting techniques and tools.

Equipment such as X-Ray machine, CT Scanners, Ultrasound scanner, Endoscope, Automatic Laboratory Analyzer etc. need specialized engineers and technicians who have been trained to maintain this specific equipment. Additionally, specialized tool/equipment may be required for testing and calibration of such equipment. This calibre of Technician/Engineer is normally employed by third party companies or vendors representing the manufacturer of the equipment.

This manual focuses on level 1 (user or first-line maintenance) and level 2 (engineer/technician level) maintenance. For purposes of this manual, level 3 maintenance may be provided by the RW Engineers/Technicians when they have the competence except for cases where Maintenance Contracts exist.

7.6 Medical Equipment that can be maintained by RWs or by outsourcing the service

The table below shows the minimum list of equipment to be maintained by the RWs and scope of maintenance for commonly used equipment in health facilities in Uganda:

S/N	Equipment Name	Level of Maintenance handled by RWs¹	Level of outsourced Maintenance from Service Providers
1	Air Compressor	Level 2	Level 3
2	Air Conditioner	Level 2	Level 3
3	Anaesthesia Unit	Level 2	Level 3
4	Analyzer, Biochemistry	Level 2	Level 3
5	Analyzer, CD4	Level 2	Level 3
6	Analyzer, Haematology	Level 2	Level 3

¹ See Section 7.5 for Maintenance levels.

S/N	Equipment Name	Level of Maintenance handled by RWs ¹	Level of outsourced Maintenance from Service Providers
7	Analyzer, Microbiology	Level 2	Level 3
8	Analyzer, Sputum	Level 2	Level 3
9	Analyzer, TB	Level 2	Level 3
10	Audiometer	Level 2	Level 3
11	Autoclave, Electric	Level 2 & 3	
12	Autoclave, Externally Heated	Level 2 & 3	
13	Balance, Analytical	Level 2 & 3	
14	Beam Balance	Level 2 & 3	
15	Bed, Gynaecology	Level 2 & 3	
16	Biological Safety Cabinet	Level 2 & 3	
17	Blood Gas Analyser	Level 2	Level 3
18	BP Machine, Aneroid	Level 2 & 3	
19	BP Machine, Digital	Level 2 & 3	
20	BP Machine, Mercury	Level 2 & 3	
21	Centrifuge	Level 2 & 3	
22	Colorimeter	Level 2	Level 3
23	Cupboard, Instrument	Level 2 & 3	
24	Cylinder Manifold System, Oxygen	Level 2 & 3	
25	Deep Freezer	Level 2 & 3	
26	Delivery Bed	Level 2 & 3	
27	Dental Chair	Level 2 & 3	
28	Dental Unit	Level 2	Level 3
29	Differential Counter	Level 2	Level 3
30	Disinfector Boiler	Level 2 & 3	
31	Drug Cabinet	Level 2 & 3	
32	ECG, 12 Channel	Level 2	Level 3
33	Electroencephalography (EEG)	Level 2	Level 3
34	Electrosurgical Unit/ Diathermy	Level 2	Level 3

S/N	Equipment Name	Level of Maintenance handled by RWs ¹	Level of outsourced Maintenance from Service Providers
35	Endoscopy – Gynaecology	Level 2	Level 3
36	Fetoscope, Doppler	Level 2 & 3	
37	Fire Extinguisher	Level 2 & 3	
38	Fume Cupboard	Level 2 & 3	
39	Generator, Petrol/ Diesel	Level 2	Level 3
41	GeneXpert Machine	Level 2	Level 3
42	HB Meter	Level 2 & 3	
43	Hot Air Oven	Level 2 & 3	
44	Incinerator	Level 2	Level 3
45	Infant Incubator	Level 2 & 3	
46	Laboratory Incubator	Level 2 & 3	
47	Medicine Trolley	Level 2 & 3	
48	Microscope Binocular	Level 2 & 3	
49	Microscope Binocular, Fluorescent	Level 2	Level 3
40	Mortuary Body Fridge	Level 2	Level 3
50	Operating Microscope	Level 2 & 3	
51	Operating Table, Electric/Manual	Level 2 & 3	
52	Oxygen Concentrator	Level 2	Level 3
53	Oxygen Therapy Apparatus	Level 2 & 3	
54	Patient Bed, Adult/ Paediatric	Level 2 & 3	
55	Phototherapy Unit	Level 2 & 3	
56	PIMA CD4 Machine	Level 2 & 3	
57	PSA Oxygen Plant	Level 2	Level 3
58	Pulse Oximeter	Level 2	Level 3
59	Radiant Infant Warmer	Level 2 & 3	
60	Refrigerator, Gas/ Electric/Kerosene	Level 2 & 3	
61	Resuscitator Adult/Infant and Crash Cart	Level 2 & 3	

S/N	Equipment Name	Level of Maintenance handled by RWs ¹	Level of outsourced Maintenance from Service Providers
62	Roller Mixer	Level 2 & 3	
63	Slit Lamp	Level 2	Level 3
64	Solar Refrigerator, Vaccine	Level 2 & 3	
65	Solar System	Level 2 & 3	
66	Spectrophotometer	Level 2	Level 3
67	Suction Machine, Electric/Manual	Level 2 & 3	
68	Syringe Pump	Level 2	Level 3
69	Ultrasonic Nebulizer	Level 2 & 3	
70	Ultrasound Scanner	Level 2	Level 3
71	Universal Anaesthesia Machine	Level 2	Level 3
72	Water filter	Level 2 & 3	
73	Water pump	Level 2 & 3	
74	Water Softener	Level 2	Level 3
75	Weighing Scale, Adult/Paediatric	Level 2 & 3	
76	Wheel Chair, Adult/Paediatric	Level 2 & 3	
77	X-ray Basic Radiography system	Level 2	Level 3
78	X-ray Film Drier	Level 2 & 3	
79	X-ray Unit, C-Arm	Level 2	Level 3
80	X-ray Unit, Fluoroscopy	Level 2	Level 3

7.7 Setting Up a Maintenance System

All RWs will strive to implement PPM. This will require an up to date medical equipment inventory at all times.

Each RW shall maintain a computerized medical equipment inventory database for each health facility using NOMAD or any other approved software by MoH. At the minimum, each RW shall

maintain a Microsoft Excel equipment inventory list extracted from the NOMAD database using the “inventory list report” tool. This list shall be extracted and maintained every quarter.

All relevant information about the equipment must be entered in the NOMAD database, including its location, serial/identification number, manufacturer details and maintenance record. A reference number should be engraved on each equipment. Additional information that should be captured includes:

- Whether equipment is maintained in-house, or
- Maintained by external agency or manufacturer

7.7.1 Routine Maintenance Procedures

The PPM activities on equipment shall be in accordance to the SOPs where available, or against a check list developed by the RW maintenance team. The specific work to be carried out may be based on guidance from the manufacturer’s maintenance/service manual.

- The frequency of maintenance shall be based on the manufacturer’s recommendation or usage.
- A heavily used equipment must be cleaned and checked more frequently than one which is used less often. Minimum standards must be met to ensure safety and reliability.

7.7.2 Engineering Personnel

Maintenance of equipment will be undertaken by trained Engineers/Technicians. In accordance with the approved MoH structure, the following maintenance personnel shall constitute the maintenance team in the RRH/RW:

- Engineer (Biomedical, Electrical, Mechanical)

- Assistant Engineering Officer (Biomedical, Electrical, Mechanical, Civil)
- Engineering Technician (Biomedical, Electrical, Mechanical)
- Engineering Assistant (Electrical, Mechanical & Civil)
- Artisan (Electrical, Plumber, Mechanical, Carpenter)

7.7.3 Reminder System

A reminder system will be put in place, so that staff are prompted to carry out tasks when they are due. A card index /sticker/calendar system or a computer programme may be used. Where card/sticker system is adopted, it will be placed on the equipment in such a way that it is visible. – See Annex-11 for Maintenance Service Sticker Form

7.7.4 Surveillance

After the programme has been set up, each RW will put in place a periodic surveillance system to ensure that records are legible and that all entries are being made. Copies of Job cards and index cards shall be stored near the equipment.

7.7.5 Standard Maintenance Work Formats

Maintenance records shall be collected and maintained using standardized formats that facilitate easy compilation and analysis of information. The format will include the following information;

- 1) Reference ID number
- 2) Equipment Name
- 3) Manufacturer
- 4) Serial Number
- 5) Date of installation
- 6) Date of maintenance
- 7) Date for next planned maintenance

- 8) Spare parts used for carrying out maintenance
- 9) Equipment functional condition before and after carrying out maintenance
- 10) Remarks on functional status

See Annex 8 for the Job Card Form

7.7.6 Special Test Equipment

Each RW will have a range of test equipment and tools to check the correct functioning of medical equipment and its compliance with electrical and other safety standards.

All RWs shall develop capacity for equipment electrical safety testing and calibration of the following parameters; (i) mass (up to 2000g), (ii) temperature, (iii) rotation/vibration speed, (iv) oxygen concentration (%) and (v) pressure.

7.7.7 Technical Library

Each RW shall maintain a fully stocked technical library with manufacturers' maintenance manuals for various medical equipment and other relevant Clinical/Biomedical engineering literature and publications.

Chapter 8: Periodic Maintenance Checklist for Common Equipment (31 types of equipment)

1) Sphygmomanometer [Blood Pressure (BP) Machine]

Function

Blood pressure is an indicator of several diseases as well as of general health. It is an easy screening test using the BP machine. A sphygmomanometer can be used to measure the blood pressure at the high point (systolic) and low point (diastolic) of the cardiac pressure cycle. Pressure is usually measured using a cuff on the upper arm.

How it works

The cuff on the arm is inflated until blood flow in the artery is blocked. As the cuff pressure is decreased slowly, the sounds of blood flow starting again can be detected. The cuff pressure at this point marks the high (systolic) pressure of the cycle. When flow is unobstructed and returns to normal, the sounds of blood flow disappear. The cuff pressure at this point marks the low (diastolic) pressure.

Pressure can be measured using a meter with dial (aneroid type), a mercury column or an electronic display. The sounds are normally detected using a stethoscope, but electronic equipment use an automatic technique with pressure sensors. The two methods do not always give the same results and the stethoscope method is generally more accurate for all types of patient.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none"> • Check whether equipment is safely packed • If mercury is spilled, seal unit and send to the technician • Remove all dust and dirt with a damp cloth
Visual checks	<ul style="list-style-type: none"> • Ensure all parts are present and are tightly fitted • Check that the display is zero when the cuff is deflated • Remove or replace any cracked rubber parts
Function Checks	<ul style="list-style-type: none"> • Before use, check that the pressure rises and returns to zero • Check correct operation of inflation bulb and valves • Remove any batteries if not in use for more than one month • Inflate to 200 mmHg and check leakage is not faster than 2 mmHg in 10 seconds

Troubleshooting

	Fault	Possible causes	Solution
1	Mercury leakage or mercury NOT at zero level	Mercury leakage or overfilling	Correction to be done by a technician
2	Mercury is dirty	Oxidation of mercury	Cleaning to be carried out by a technician
3	Pressure does NOT increase easily or Pressure increases after inflation	Valve or tube blockage	Remove and clean all valves and tubes. Reassemble and test
4	Aneroid instrument does NOT return to zero	Zero setting has moved	Rotate collar on the base until zero setting is achieved and tighten.
5	Pressure does NOT remain steady	Leakage of air	Isolate leak by closing off parts of tubing. Replace leaking section and retest

2) Stethoscope

Function

A stethoscope is used to listen to sounds within the body. These might be sounds generated by breathing, coughing, blood flow or the stomach. The sounds are picked up and transmitted to the ears of the medical staff for diagnosis.

How it works

A membrane on the stethoscope head picks up the vibrations caused by internal sounds and transmits them to the stethoscope tube. The sounds pass up the tube through the earpiece to the user. The stethoscope head also contains an open bell which is used to pick up lower frequency sounds. The head picks up the sound from a wide area so it sounds loud to the user. Care must therefore be taken not to hit or shout into the stethoscope while in use.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none"> • Check that the equipment is safely packed • Remove any visible dirt • Remove all dust and dirt with a damp cloth Remove earpieces and clean the inside with warm water
Visual checks	<ul style="list-style-type: none"> • Ensure all parts are present and are tightly fitted • Remove or replace any cracked rubber parts • Replace membrane if broken
Function Checks	<ul style="list-style-type: none"> • Tap gently before use to check operation • Check tube holder rotates easily within headpiece • Check sound can be heard from both sides of headpiece

Troubleshooting

	Fault	Possible causes	Solution
1	Faint or NO sound heard	Leakage or blockage	Remove all parts and check for leakage and blockage.
2	Tube connector does NOT stay in headpiece	Broken locking mechanism	Repair to be done by a technician
3	Parts damaged or faulty	Broken part	Replace with part taken from other units

3) Refrigerator (Electric and Solar)

Function

A refrigerator is a hospital plant that keeps health facility vaccines, reagents and blood at a required environment (e.g. temperature and humidity).

How it works

An electrical refrigerator can be a compression or absorption type. Compression fridges use electrical power as a source of power and it has four major components;

- 1– A compressor,
- 2– A condenser,
- 3– An Expansion valve/capillary tube,
- 4– An Evaporator

A compressor has a pump and electrical coil, both inside the housing. The coil gets electrical power to form an electrical field to drive the pump which then pumps the refrigerant in form of gas with a high pressure to the condenser.

The condenser receives the refrigerant, removes heat and condenses refrigerant into a liquid. The liquid refrigerant still at high pressure goes to an expansion valve/capillary tube.

The Expansion valve or capillary tube drops the refrigerant pressure and temperature and this is where cooling starts.

The Evaporator is the cooling part of the refrigerator which evaporates the liquid refrigerant back into gas form and removes heat from the room.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none"> • Use a damp cloth to clean the solar panel and fridge body. • Use a damp cloth to clean evaporator after melting ice. Never use a sharp object to remove ice. • Use a brush or dry cloth to clean the condenser. • Use a brush to clean the compressor. • Use a brush to clean the charge regulator.
Visual checks	<ul style="list-style-type: none"> • Check door rubbers if not damaged. • Check thermometer position in the fridge. • Check room ventilation where the fridge is located. • Check discharge level on the charge regulator. • Check ice formation inside the fridge. • Check dust on solar panels. • Check if solar panels are not shaded.
Function Checks	<ul style="list-style-type: none"> • Check whether the power source is available. • Check thermostat position inside the fridge. • Check thermometer reading from the fridge. • Check fridge stability on the ground.

Troubleshooting

	Fault	Possible causes	Solution
1	Compressor NOT running.	No power supply. Burnt compressor.	<ul style="list-style-type: none"> • Check power supply. • Replace the compressor.
2	Compressor running but NOT cooling.	Gas leakage. Weak compressor pump. Chocked drier.	<ul style="list-style-type: none"> • Braze pipe, charge the fridge. • Replace compressor. • Replace filter drier.
3	Fridge over freezing.	Faulty thermostat. Wrong Thermostat setting.	<ul style="list-style-type: none"> • Replace thermostat. • Set thermostat properly.
4	Too much ice formation inside fridge.	Evaporator fan not running.	<ul style="list-style-type: none"> • Replace evaporator fan.

4) Refrigerator (Gas and Kerosene)

Function

The function of a refrigerator is to bring down the room temperature to the required one.

How it works

The absorption system is different from the compression system. It uses heat energy instead of mechanical energy to make a change in the conditions necessary to complete a refrigeration cycle. The system may use butane, kerosene, steam or electrical heating element as a source of heat and ammonia solution as a refrigerant.

The system has few moving parts, small systems have moving parts only in the heat source valves and controls which are used. Larger systems also use circulating pumps and fans as moving parts.

Absorption systems also have 4 major components. 1–Generator, 2–Condenser, 3–Absorber/Receiver and 4–Evaporator. In this case the compressor is replaced by a heater and Generator.

Absorption system uses a generator charged with water and ammonia, a heat source usually kerosene flame heats this solution in the generator, the ammonia becomes vaporized and is driven off to the condenser, and the condenser at the top of the system condenses the ammonia vapor into liquid. The liquid flows by gravity into the liquid receiver and then into the evaporator. During the generating cycle, little or no refrigerating effect is taking place. As the system cools, pressure drops causing the liquid ammonia in the evaporator to boil (cooling) and absorb heat from the room, the cycle is complete when vaporized ammonia is absorbed back to the generator.

Maintenance Checklist

Cleaning	<ul style="list-style-type: none"> • Use a damp cloth to clean the fridge body. • Use a brush to clean the condenser. • Use a brush to clean the generator. • Use a damp cloth to clean the inside of the fridge when it is off and ice is melting.
Visual checks	<ul style="list-style-type: none"> • Check if gas cylinder is open and there is no leakage. • Check if spirit level is in a better position/leveling. • Check if thermostat probe is in a better position.
Function Checks	<ul style="list-style-type: none"> • Check if the burner or nozzles are clean and okay. • Check if gas cylinder is open and has gas. • Check if thermostat is functioning well. • Check the level of ice formation on evaporator.

Troubleshooting

	Fault	Possible causes	Solution
1	NOT cooling well.	<ul style="list-style-type: none"> • Nozzles blocked. • Too much ice formation. 	<ul style="list-style-type: none"> • Clean nozzles. • Carry out manual defrosting.
2	Over cooling.	<ul style="list-style-type: none"> • Thermostat not cutting off gas or heating element. 	<ul style="list-style-type: none"> • Replace the thermostat.
3	Burning Gas available but NOT cooling.	<ul style="list-style-type: none"> • Ammonia leakage. 	<ul style="list-style-type: none"> • Braze the pipe and charge the system.
4	Generator heating okay, No leakage but no cooling.	<ul style="list-style-type: none"> • Bubble formation in the system. 	<ul style="list-style-type: none"> • Turn the system horizontal for at least two days.

5) Weighing Scales

Function

Measuring patient weight is an important part of monitoring health as well as calculating drug and radiation doses. It is therefore vital that scales continue to operate accurately. They can be used for all ages of patient and therefore vary in the range of weights that are measured. They can be arranged for patients to stand on, or can be set up for weighing wheelchair bound patients. For infants, the patient can be suspended in a sling below the scale or placed in a weighing cot on top of the scale.

How it works

Mechanical scales have a spring deflected by patient weight. The spring pushes a pointer along a display or rotates a disc to indicate weight. Electronic scales have a sensor (pressure sensor) that bends under patient weight and the circuitry converts this to displayed digits.

This pressure sensor under strain induces an electric potential which is directly proportional to the applied force (weight).



Baby weighing scale, dial type



Baby weighing scale, spring type



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none"> • Wipe off dust and replace dust cover after checks • Clear away any dirt or hair on controls and feet • Clean exterior with damp cloth and dry off • Clean off then repaint any exposed or rusted metal
Visual checks	<ul style="list-style-type: none"> • If bent, cracked or damaged, send for repair • Tighten any loose screws and check parts are fitted tightly
Function Checks	<ul style="list-style-type: none"> • Check zero at start of day and before each patient • Check reading is accurate using a known weight • Send for repair if inaccurate or sticking • Replace battery if display shows low battery

Troubleshooting

	Fault	Possible causes	Solution
1	Zero point can NOT be set	Scale is not level Zero control broken or internal part jammed	Set scale on level ground and Retest Repair or change parts
2	Movement is stiff or jerky	Dirt lodged inside Internal blockage	Remove any visible dirt or foreign body and reset Lubricate /repair
3	Reading is inaccurate	Zero not properly set Calibration error	Reset to zero Recalibrate
4	Electronic display is blank	Battery / power failed	Replace battery or power supply and retest

6) Suction Machines (Aspirators)

Function

Suction machines (also known as aspirators) are used to remove unwanted fluid from body cavities. They are found in operating theatres, delivery suites, ENT and emergency departments. Smaller specialized suction machines are used in dental departments.

How it works

Suction is generated by a pump. This is normally an electrically powered motor, but manually powered versions are also often found. The pump generates a suction that draws air from a bottle. The reduced pressure in this bottle then draws the fluid from the patient via a tube. The fluid remains in the bottle until disposal is possible. A valve prevents fluid from passing into the motor itself.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none">• Wipe dust off exterior of the equipment and bottle cover.• Wash bottle and patient tubing with sterilizing solution• Wipe round bottle seal with damp cloth, replace if cracked• Remove dirt from wheels / moving parts
Visual Checks	<ul style="list-style-type: none">• Check if all fittings & accessories are mounted correctly• Check if filter is clean• Check if parts are fitted and replace any cracked tubes• Check if mains cable has no bare wire and is not damaged
Function Checks	<ul style="list-style-type: none">• Check that all switches and vacuum control operate correctly• Check for air and liquid leakage

Troubleshooting

	Fault	Possible causes	Solution
1	Machine is NOT running	No power from mains socket Fuse blown Electrical cable fault	Check power switch is on. Check mains power. Check for leaks or short circuits causing the fuse to blow and correct this. Replace fuse. Try cable on another piece of equipment to rule out internal cable damage
2	Poor fluid flow, pressure gauge low	Tube /seal / bottle leaking or Disconnected Air outlet valve blocked Control valve stuck Internal or control error	Close different tubes by bending. When pressure gauge changes, leakage point has been passed. Replaced damaged tube or seal. Clean outlet valve Operate control valve through full range. Send for repair if it is stuck. Refer to a technician
3	Poor fluid flow, pressure gauge high	Blocked filter or tube	Disconnect each tube one at a time. When air flow is stopped, blockage has been passed. Replace filter or unblock tube.
4	Filter discolored	Floating valve broken	Change filter, clean or replace floating valve
5	Electrical shocks	Wiring fault	Refer to an electrician
6	Manual suction is jammed	Internal slider stuck	Refer to a technician for greasing

7) Ultrasonic Nebulizers

Function

Is a device used to administer medication in the form of a mist inhaled into the lungs. Nebulizers are commonly used for treatment of asthma and other respiratory diseases. The reason for using a nebulizer for medicine to be administered directly to the lungs is that small aerosol droplets can penetrate into the narrow branches of the lower airways. Large droplets would be absorbed by the mouth cavity, where the clinical effect would be low.

How it works

Ultrasonic Nebulizers use ultrasonic power as means to break up medical solutions or suspensions into small droplets, these small droplets are passed for direct inhalation either through the mouthpiece of the device or a hose set. An Ultrasonic Nebulizer uses a small crystal to generate vibrations in the solution that cause droplets to break off.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none">• Clean and sterilize mouthpiece and medicine chamber• Wipe dust from the machine and replace cover after checks
Visual Checks	<ul style="list-style-type: none">• Check all parts are present and tightly fitted• Check all moving parts move freely, all holes are unblocked• If mains plug, cable or socket are damaged, replace them• If chamber and tube seals are damaged, replace them.

Function Checks	<ul style="list-style-type: none"> • Check the whole system functions before use • Before next use, check that there is adequate nebulization. • Check that the compressor fan is working without excessive noise.
-----------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Troubleshooting

	Fault	Possible causes	Solution
1	Equipment is NOT working	No power from mains socket/ Blown mains fuse Electrical cable fault	Check power switch is on. Replace fuse with correct current ratings. Check mains power. Try cable on another piece of equipment.
2	Equipment is on but flow is absent	Filter is blocked Pipe is twisted or nebulizer chamber/mouth piece is blocked.	Clean filter Connect pipe properly, clean chamber / mouthpiece
3	Inadequate nebulizing amount	Output adjustment not correctly set. Mouth piece cracked Vibration generator weak	Adjust output as directed in the user manual Replace mouthpiece Replace vibration generator.
4	Electrical shocks	Improper earthing/grounding of the machine. Bare wires touching the body of the machine	Earth the machine appropriately Insulate all live conductors


8) Oxygen Concentrators

Function

An oxygen concentrator draws in room air, separates the oxygen from the other gases in the air and delivers the concentrated oxygen to the patient. When set at a rate of two litres per minute, the gas that is delivered by the concentrator is more than 90% oxygen. It is used for situations where bottled gas supply is impractical or expensive, and can be used by patients in the hospital or the home.

How it works

Atmospheric air consists of approximately 80% nitrogen and 20% oxygen. An oxygen concentrator uses air as a source of oxygen by separating these two components. It utilizes the property of zeolite granules to selectively absorb nitrogen from compressed air. Atmospheric air is gathered, filtered and raised to a pressure of 20 psi (138kPa) by a compressor. The compressed air is then introduced into one of the canisters containing zeolite granules where nitrogen is selectively absorbed leaving the residual oxygen available for patient use. After about 20 seconds the supply of compressed air is automatically diverted to the second canister where the process is repeated enabling the output of oxygen to continue uninterrupted. While the pressure in the second canister is at 20 psi (138kPa) the pressure in the first canister is reduced to zero. This allows nitrogen to be released from the zeolite and returned into the atmosphere. The zeolite is then regenerated and ready for the next cycle. By alternating the pressure between the two canisters, a constant supply of oxygen is produced and the zeolite is continually being regenerated. Individual units have an output of up to five litres per minute with an oxygen concentration of up to 95%.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none"> Remove any dust / dirt with a damp cloth and dry off Fill humidifier bottle up to the level marker with clean distilled water Wash filter in warm water and dry. Replace if damaged Clean humidifier bottle thoroughly and dry off
Visual Checks	<ul style="list-style-type: none"> Check all screws, connectors, tubes and parts tightly fitted Replace humidifier bottle if covered with limescale. Check mains plug, cable or socket for any damage and replace if necessary.

Function Checks	<ul style="list-style-type: none"> • Check oxygen flow before setting up line to a patient. • Run machine for two min. and check that no alarm goes off. • Check functionality of flow meter. Oxygen flow rate should vary with flow control (see bubbles).
-----------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Troubleshooting

	Fault	Possible causes	Solution
1	Unit NOT operating, power failure alarm sounds	No power from mains socket Concentrator circuit breaker has been set off.	Check mains switch is on and cable inserted. Replace fuse. Press reset button if present
2	Unit NOT operating, no power failure alarm	Alarm battery dead	Replace battery and test as above
3	No oxygen flow	Flow not visible Tubes not connected tightly Water or dirt blocking the oxygen tubing Blocked flow meter or humidifier bottle	Place tube under water and look for bubbles. If bubbles emerge steadily, gas is indeed flowing Check tubing and connectors are fitted tightly Remove tubing, flush through and dry out before replacing Replace meter / bottle or refer to a technician
4	Temperature light or low oxygen alarm is on	Unit over heated or obstructed	Remove any obstruction caused by drapes, bedspread, wall Clean filters. Turn unit off, using standby oxygen system. Restart unit after 30 min.
5	Electric shocks	Wiring fault	Refer to an electrical technician
6	Concentrator is humming but can NOT start	Defective capacitor	Test with a multi-meter and replace
7	Service warning light persists for more than 5 min with or without audible alarm	Dirty filters	Clean or replace the filters

9) Oxygen Cylinders and Flow meters

Function

Medical gases such as oxygen, nitrous oxide etc. are intended for administration to a patient in anaesthesia, therapy or diagnosis. An oxygen cylinder is a cylindrically shaped metal container used to store oxygen that has been compressed to a very high pressure. Oxygen cylinders, which come in different sizes, are usually coloured black with a white top; in some cases, it may be a small cylinder that is entirely black. The black colour helps to differentiate it from other gases that are stored in similar containers. Cylinders are fitted with customized valves (either bull nose or pin index type) with valve guards, which are opened with valve keys.

A flow meter is an instrument used to measure the flow rate of a liquid or a gas. In healthcare facilities, gas flow meters are used to deliver oxygen at a controlled rate either directly to patients or through medical devices. Oxygen flow meters are used on oxygen tanks and oxygen concentrators to measure the amount of oxygen reaching the patient or user. Sometimes bottles are fitted to humidify the oxygen by bubbling it through water.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none">• Ensure delivery tubes and masks are sterile• If humidifier bottle is used, refill with clean water• Clean cylinder, valve and flow meter with a damp cloth
Visual Checks	<ul style="list-style-type: none">• Check that the cylinder is the correct type and marked oxygen.• Check all parts are fitted tightly and correctly• Check for leakage: hissing sound or reduction in pressure

Function Checks	<ul style="list-style-type: none"> • Before use, ensure cylinder is filled & flow is present • Close cylinder valve after each use. • Remove valve dust with brief, fast oxygen flow • Check flow can be varied using flow control
-----------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Troubleshooting

	Fault	Possible causes	Solution
1	No oxygen is flowing	Cylinder is empty Flow meter knob or cylinder valve is closed. Faulty regulator	Replace cylinder Open valves, then check flow meter registers flow Close all valves and replace regulator
2	Leakage from cylinder or Flow meter	Cylinder is not connected to regulator properly Faulty or missing washer Flow meter seal damaged or loose Cylinder faulty	Tighten all fittings Replace washer Tighten flow meter Label "Faulty" and return to manufacturer
3	Leakage can NOT be located	Leakage too small to be heard	Apply detergent solution (NOT oily soap) to joints. Bubbles will show at leak point. Clean/replace washer and tighten at that joint.
4	Flow meter ball NOT moving, yet oxygen is flowing	Faulty flow meter	Close all valves, disconnect flow meter and clean inside. Reconnect and test.
5	Pressure gauge does NOT show pressure, yet oxygen is flowing	Faulty pressure gauge	Refer to biomedical technician for Replacement

10) General X-Ray Machine

Function


X-Ray machines are used for imaging bones and hard tissues and diagnosing fractures, joint defects, choked lungs etc. Sometimes contrast agents are also used to highlight any defects in the abdomen under X-rays.



How it works

X-rays are high energy electromagnetic waves. The transformer produces a high voltage that directs electrons onto a target in the machine head. X-rays are produced by the target and are directed into beams by a collimator towards the human body. Soft body tissue absorbs less X-rays, i.e., more radiation goes through, whereas bone and other solids prevent most of the X-rays from going through. A photographic film or electronic sensor displays how much X ray has passed through, forming an image of the interior of the body. Bone appears nearly white, because few X-rays strike the corresponding part of the film, leaving it largely unexposed; soft tissue allows much more radiation to pass through, darkening the film in those places.

Users must ensure proper radiation safety protocols and supervision is in place.



(Control panel and transformer not shown)

Maintenance Checklist

Care and Handling	<ul style="list-style-type: none"> • Remove any items or foreign body from equipment • Clean all dust and dirt from the X-Ray machine and room
Visual Checks	<ul style="list-style-type: none"> • Check that all parts are present and connected • Check that cables are not twisted and remove from service if any damage is visible • If any plug, cable or socket is damaged, refer to distributor • Check lead aprons for any defects • Check table, cassette holder & grids for smooth movement
Function Checks	<ul style="list-style-type: none"> • Switch on power and check all indicators function • If machine has not been in use, wear lead apron and check whether exposure indicator lights on switch operates • Check collimator bulb, replace with correct type if needed • Check power supply voltage, and regulation • Check protection earthing for the Equipment • Check that internal earth connections are firm

Troubleshooting

	Fault	Possible causes	Solution
1	X-Ray unit does NOT switch on	Mains power not connected	Check the machine is plugged into the mains socket and that all switches are on. Replace fuse. Check if mains power is present at socket using equipment known to be working. Rewiring by electrician.
2	X-Ray machine NOT exposing, even when power is on.	Safety interlock is on Exposure switch cable problem Internal error	Check safety locks, all switches Check for any loose connection Refer to distributor's technician
3	Poor X-Ray image quality	X-Ray tube problem	Refer to distributor's technician
4	The table does NOT move.	Table motor or cable problem. Safety switch or fuse problem Control circuit problem	Check all cable connections Check relevant fuse or switch Refer to biomedical technician
5	Electrical shocks	Wiring fault	Refer to distributor's technician immediately

11) Ultrasound Scanner

Function

Diagnostic ultrasound machines are used to give images of structures within the body. Other kinds of machines (e.g. therapeutic and lithotripsy) are not dealt with herein. The diagnostic machine probes, which produce the ultrasound, come in a variety of sizes and styles, each type being produced for a particular special use. Some require a large trolley for all the parts of the unit, while the smallest come in a small box with only an audio loudspeaker for output. They may be found in cardiology, maternity, outpatients and radiology departments and will often have a printer attached for recording images. Unlike X-rays, ultrasound poses no danger to the human body.

How it works

The ultrasound probe contains a crystal that sends out bursts of high frequency vibrations that pass through gel and then through the body. Soft tissue and bone reflect echoes back to the probe, while pockets of liquid pass the ultrasound straight through. The echoes are picked up and arranged into an image displayed on a screen. The machine offers a number of processing options for the signal and image and also allows the user to measure physical features displayed on the screen. This requires the machine to incorporate a computer.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none">• Wipe dust off exterior and cover equipment after use• Remove any items or foreign body from equipment• Wipe probe with alcohol-free tissue or cloth• Remove, clean and dry external filter if present
-------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Visual Checks	<ul style="list-style-type: none"> • Check that all fittings and accessories are mounted correctly • Check that cables are not twisted and probe is safely stored • Check that mains plug screws are tight • Check that mains cable has no bare wire and is not damaged
Function Checks	<ul style="list-style-type: none"> • If machine has not been in use, run and test briefly

Troubleshooting

	Fault	Possible causes	Solution
1	Equipment is NOT running	No power from mains socket Electrical cable fault	Check power switch is on. Replace fuse. Check mains power is working. Repair or exchange the cable. Try cable on another piece of equipment.
2	Fuse keeps blowing	Power supply or cable fault	Refer to distributor's technician
3	Probe head damaged or noisy	Possible internal fault	Exchange probe.
4	Image quality poor	Gel insufficient Controls set incorrectly Mains voltage is too low	Use more ultrasound gel Check controls for correct positioning and operation (refer to user manual) Use voltage stabilizer
5	Display / computer error	Software fault	Turn machine off and restart. If problem persists, refer to distributor's technician
6	Electrical shocks	Wring problem	Rewiring/repair by a technician

12) Electrosurgical Units (ESU) / Diathermy

Function

Electrosurgery is the application of a high–frequency electric current to biological tissue as a means to cut, coagulate, desiccate, or fulgurate tissue. Its benefits include the ability to make precise cuts with limited blood loss in hospital operating rooms or in outpatient procedures. Cautery, or electrocautery, is the application of heat to tissue to achieve coagulation. Although both methods are sometimes referred to as surgical diathermy, this chapter avoids the term as it may be confused with therapeutic diathermy, which generates lower levels of heat within the body.

How it works

In electrosurgical procedures, the tissue is heated by an alternating electric current being passed through it from a probe. Electrocautery uses heat conduction from an electrically heated probe, much like a soldering iron. Electrosurgery is performed using an electrosurgical generator (also referred to as power supply or waveform generator) and a hand piece including one or several electrodes, sometimes referred to as an RF Knife, or informally by surgeons as a "Bovie knife" after the inventor. Bipolar electrosurgery has the outward and return current passing through the handpiece, whereas monopolar electrosurgery returns the current through a plate normally under the patient. Electrosurgery is commonly used in dermatological, gynaecological, cardiac, plastic, ocular, spine, ENT, orthopedic, urological, neuro– and general surgical procedures as well as certain dental procedures.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none">• Remove any dust / dirt and replace equipment cover• Unplug, clean outside with a damp cloth and dry off
Visual Checks	<ul style="list-style-type: none">• Check that all fittings and cables are properly connected• Check there are no signs of spilled liquids or cable damage• Inspect filters, clean or replace if needed.• If any plug, cable or socket is damaged, replace
Function Checks	<ul style="list-style-type: none">• Check foot / probe switch smooth operation.• Check returns plate cable disconnection alarm before use.

	<ul style="list-style-type: none"> • Check proper operation of all controls, indicators and visual displays on the unit. • If not recently used, check operation on wet soap
--	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Troubleshooting

	Fault	Possible causes	Solution
1	Equipment is NOT running	No power from mains socket Electric cable fault	Check power switch is on. Replace fuse. Check mains power. Try cable on another piece of equipment.
2	Equipment is on but shows error signal	Footswitch pedal may have been depressed as unit is turned on or front panel buttons may be stuck. Probe, patient cable or plate malfunction	Note error code and turn unit off. Check footswitch and front panel buttons. Disconnect all foot pedals. Turn on unit again. Check connections and plugs on all cables are tight.
3	Equipment is on but output is absent, weak or intermittent	Power setting is too low Malfunctioning parts Incorrect connection	Adjust power, check manual Check connection or replace parts Check correct probe/footswitch cord are well connected
4	Continuous interference with monitors	Faulty ground connection Poor filtering systems in monitoring device	Check all monitors and power connections. Use separate outlets for each device. Replace monitoring device
5	Monitor interference occurs only when electrosurgery is activated	Metal to Metal sparking Cords and cables are bundled, touching or damaged High power setting	Check all connections are tight. Remove cable cluttering, replace damaged cords Reduce power setting, use blend Mode
6	Pacemaker or internal cardiac defibrillator interference	Equipment activation is causing battery or implant malfunction	Stop procedure immediately, perform emergency care and call implant supplier before restarting procedure
7	Electrical shocks to user	Wiring fault	Refer for rewiring by a technician

13) Tables (for Operating Theatres and Delivery)

Function

Tables are required to hold the patient in a position that is comfortable for both the patient and for medical staff during procedures. They can include dedicated supports for head, arms and legs and often have movable sections to position the patient appropriately. They are made both with wheels and on static platforms and can have movements powered by electric motors, hydraulics or simply manual effort. They can be found in emergency departments, operating theatres and delivery suites.

How it works

Where the table has movement, this will be enabled by unlocking a catch or brake to allow positioning. Wheels have brakes on the rim or axle of the wheel, while locks for moving sections will normally be levers on the main table frame. Care should be taken that the user knows which lever applies to the movement required, as injury to the patient or user may otherwise result. The table will be set at the correct height for patient transfer from a trolley then adjusted for best access for the procedure.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none">• Clean, dry and disinfect table/ underneath/ base• Wipe off any escaped oil or grease from joints
Visual Checks	<ul style="list-style-type: none">• Check that all parts are present and tightly fitted• Check that no oil is leaking• Inspect mattress and table for signs of wear• Replace any worn or damaged items

Function Checks	<ul style="list-style-type: none"> • Check essential movements before use • Check wheel brakes function and wheels rotate • Ensure all moving parts can move, applying grease if needed
-----------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Troubleshooting

	Fault	Possible causes	Solution
1	Table can NOT be relocated	<p>Wheels jammed</p> <p>Electric motor not operational (electrically driven table)</p>	<p>Clean wheels, remove obstruction</p> <p>Check power to table Replace fuse if blown If problem persists, refer to technician</p>
2	Table section or body can NOT be moved	<p>Lock or lever is jammed</p> <p>No power to electric table</p> <p>No oil in hydraulic table</p>	<p>Clean jammed part, remove rust and dirt, lightly oil and replace</p> <p>Check correct switch is used</p> <p>Check power and fuses</p> <p>Refill hydraulic oil if needed</p> <p>Check no leakage occurs</p>
3	Oil leakage from hydraulic table	<p>Oil leakage</p> <p>Worn out gasket oil seals.</p>	<p>Locate leak and block it.</p> <p>Clear spillage.</p> <p>Refer to technician.</p>
4	Electric shocks	Wiring fault	Refer to technician immediately

14) Lamps

Function

There are many kinds of sources of light used in medicine. This section deals with large lights for operating theatres or delivery suites, ultraviolet or infrared phototherapy units, ophthalmic slit lamps, handheld and head worn lamps for ENT clinics and domestic torches. However, the principles here will help in the maintenance of any kind of light source.

How it works

Each type of lamp will have a power source with switch and a bulb. Some will also have controls for the brightness or focus of the light, while others will also have lenses to direct the light where required. Some lights operate off mains electricity, while others use batteries instead. Some lights have both, using the batteries for back-up power in case of mains supply failure. Electric bulbs and batteries have limited life and will need regular checking. A stock of spares should be kept of all the correct voltages and wattages (ratings) of parts.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none">• Unplug, clean outside of the lamp r with a damp cloth and dry off.• Clean and sterilize the lamp head handle
Visual Checks	<ul style="list-style-type: none">• Check that all lamps are functioning.• Check that there are no cracks in glass / covers or liquid spillages• Tighten any loose screws and check parts are fitted

Function Checks	<ul style="list-style-type: none"> • Check that switches and focus knob operate correctly • Check working condition of moving parts
-----------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------

Troubleshooting

	Fault	Possible causes	Solution
1	No light or power on visible	No power at mains socket Dead battery Blown bulb Battery leakage Electrical cable fault Internal wiring fault	Check power switch is on. Replace fuse. Check mains power. Charge or replace battery Replace bulb with correct voltage and wattage Remove batteries, clean battery terminals and replace with new battery Try cable on another piece of equipment. Refer to an Electrician
2	Fuse / bulb keeps blowing	Fuse or bulb is wrong rating Power supply or cable fault	Replace with correct rating Refer to an electrician
3	Light can NOT be made bright enough	Dirt on lens or tube Poor power supply Wrong bulb rating Control malfunction	Clean area with dry, clean cotton Check power line or replace Batteries Check bulb rating is correct Refer to an electrician
4	Electrical shocks	Wiring fault	Refer to an electrician

15) Anesthesia Machine

Function


The anesthesia machine is used to support the administration of anesthesia. The most common type of anesthesia machine is the continuous-flow type, which is designed to provide an accurate and continuous supply of medical gases (such as oxygen and nitrous oxide), mixed with an accurate concentration of anesthetic vapour (such as halothane or isoflurane), and deliver this to the patient at a safe pressure and flow. Modern machines incorporate ventilator, suction unit and patient monitoring devices.

How it works

Oxygen (O_2), nitrous oxide (N_2O) and sometimes air sources are connected to the machine. Through gas flowmeters (or rotameters), a controlled mixture of these gases along with anesthetic vapour passes through a vaporizer and is delivered to the patient. Sometimes a ventilator is also connected with the machine for rebreathing thus making it a closed circuit. With ventilators or a re-breathing patient circuit, soda lime canisters are used to absorb the exhaled carbon dioxide and fresh gases are added to the circuit for reuse. Pressure gauges are installed on the anesthesia machine to monitor gas pressure. Generally, 25% (or 21%) oxygen is always kept in the circuit (delivered to patient) as a safety feature. The device which ensures this minimum oxygen in the circuit is called a hypoxic guard. Some basic machines do not have this feature,

but have a nitrous lock which stops the delivery of N_2O in absence of O_2 pressure.

Machines give various alarms to alert operators.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none"> Remove water and waste matter from inside Clean inside and outside with a damp cloth and dry off
Audio–Visual checks	<ul style="list-style-type: none"> If any leak is audible, check with soapy solution Check all seals, connectors, adapters and parts are tight Check all moving parts move freely, all holes are unblocked Check connections for leakage with soap solution & dry off Check all fittings for proper assembly Replace soda lime if it has turned blue Replace any deteriorated hoses and tubing If seal, plug, cable or socket are damaged, replace
Function Checks	<ul style="list-style-type: none"> After use, depressurize system and replace all caps / covers Before use, check pressure gauges rise and there are no leaks

Troubleshooting

	Fault	Possible causes	Solution
1	Equipment NOT running	No power at mains socket Electric cable fault	Check power switch is on. Replace fuse. Check mains power is present at socket using equipment known to be working. Rewiring if power not present.
2	No Gas output	No O ₂ pressure in cylinder / gas supply Check pressure gauges (4 bar or 4 kg/cm ²)	Restore gas supply or replace gas cylinders Replace O ₂ cylinder and/or N ₂ O cylinder in case of low pressure.
3	O ₂ failure alarm NOT working	Alarm battery is low. Alarm device is not working	Check alarm setting range, alarm on/off, alarm lamp blow out
4	Machine has leaks	Poor seal (commonly occurring around tubing connections, flow valves and O ₂ /N ₂ O yokes) Cylinders not seated properly	Clean leaking seal or gasket, replace if broken. Refit cylinders in yokes and retest. If leaks remain, repair by technician.
5	Flow meter fault	Over tightening of the needle valve or sticking of the float / ball	Repair by technician or outsourcing service
6	Electrical shocks	Wiring fault	Rewiring/repair by electrician immediately

16) Infant Incubator

Function

Infant incubators are classified into closed type and open type. An infant incubator is intended to hold a neonate born prematurely or with some disabilities for life support. The closed type of incubator is applied in the case of treatment of neonate jaundice. The open type is often used to facilitate the surgical treatment for neonates who have undergone surgical operation. Humans feature excellent homeostasis of body temperatures. Since homeostasis of neonates is lower, keeping of the body temperature is essential for the neonates. Premature and impaired neonates have further lower homeostasis and thus incubators are indispensable for sustaining their life.

How it works

In practice, incubators cannot be applied in the environment where the outside air temperature exceeds the body temperature. To keep a certain temperature, every incubator is furnished with a heater and a fan motor for air circulation. For the motor, hour meter management is essential.

Note: The hood in the closed type incubators is often made of acrylic material, therefore alcohol and strong chemical agents should not be applied to the hood.



Maintenance Checklist

Cleaning	<ul style="list-style-type: none">• Unplug, clean outside with damp cloth and dry off• Clean the access ports and temperature sensor• Wash/replace the air filters, dry thoroughly for reuse• Disinfect water chamber of humidifier
Visual Checks	<ul style="list-style-type: none">• Check all fittings and accessories are mounted correctly• Check mains cable has no bare wire and is not damaged• Check doors, cable and tray. Repair if damaged

Function Checks	<ul style="list-style-type: none"> • Drain off the water tray. Run machine for 30 minutes to dry the tray. Refill tray with sterile water before use. • Check all controls operate correctly • Check the readings of thermometer and oxygen sensors change when breathed upon • Check that batteries are working properly.
-----------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Troubleshooting

	Fault	Possible causes	Solution
1	Incubator is NOT running	No power from mains socket Electrical cable fault	Check power switch is on. Replace fuse. Check mains power. Try cable on another piece of equipment. Repair by electrician, if required.
2	Fuse keeps blowing	Power supply or cable fault	Check the Manual and troubleshoot properly
3	Alarms NOT working	Alarm battery dead. Alarm off setting	Replace the battery and recheck. Check alarm setting.
4	Temperature NOT properly controlled	Temperature probe and sensor not working Incubator placed indirect sunlight or near a fan. Fan or air duct problem	Check the temperature probes and sensor connections. Replace the temperature probe or sensor, if required. Move incubator if placed near heat or draught Unblock air duct if obstructed.
5	Incubator NOT heating even when the heater lamp is on.	Heating element problem	If accessible, replace heating element.

17) Infant Warmer

Function

An infant warmer system provides a controlled source of radiant heat for infants and pediatric patients.

How it works

The heater assembly consists of a radiant heater controlled by electronic circuits. A probe is used to monitor the patient's skin temperature. The patient's skin temperature is continuously displayed. Alarms activate to alert the operator of a low or high patient temperature, a skin temperature probe failure, a power failure, equipment failure or a check patient prompt.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none">• Clean the warmer at least once a week or after each patient.• The mattress, X-ray tray, bed and side panels may be cleaned without immersing by using a disinfecting agent safe for use on the materials.
Visual Checks	<ul style="list-style-type: none">• For units with casters, check that all casters are in firm contact with the floor and that the warmer is stable and moves freely.• Examine the unit for objects placed on top of the heater assembly.
Function Checks	<ul style="list-style-type: none">• Check that all accessories are mounted securely and that the load limits are not exceeded.• Connect the warmer to power source and verify the operation of the control panel.• Check operation of phototherapy warning light.

Troubleshooting


	Fault	Possible causes	Solution
1	Alarm sounding and manual indicator flashing.	Software has detected a fault	Use information mode to obtain error code.
2	Manual indicator flashing accompanied by monotone alarm.	Thermal cut out in head has switched open circuit due to overheating.	Reset thermal cut– out and identify cause.
3	Manual indicator flashing accompanied by rapid pulsing two tone alarm.	Fault in the control or power board.	Replace defective PCB.
4	Power failure	Power switch is off, internal fuses, mains lead and internal harnesses defective.	Switch on power, replace fuse, defective part or PCB.
5	Wrong Skin temperature reading	Skin sensor is damaged, not inserted fully or not the right type.	Replace skin sensor.
6	Unable to provide stable control of skin temperature when in baby/patient mode	<p>Skin sensor is poorly attached to patient or heat path between patient and element is disrupted.</p> <p>Alarm is sounding and mute button is not being pressed.</p>	<p>Attach skin sensor in the correct position and ensure there is no heat path disturbance.</p> <p>Press mute when warmer alarms. This ensures heat is restored to the patient.</p>

18) Autoclaves

Function

Sterilization is the killing of micro-organisms that could harm patients. It can be done by flame or (steam, air, flame or boiling) or by chemical means. Autoclaves use high pressure steam and Sterilizers use boiling water mixed with chemicals to achieve this. Materials are placed inside the unit for a carefully specified length of time. Autoclaves achieve better sterilization than boiling water sterilizers.

Heat is delivered to water either by electricity or flame. This generates high temperature within the chamber. The autoclave also contains high pressure when in use, hence the need for pressure control valves and safety valves. Users must be careful to check how long items need to be kept at the temperature reached.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none"> • Remove water and waste matter from inside • Clean outside of equipment and power cord • Drain water reservoir and clean inside of reservoir using vinegar solution
Visual Checks	<ul style="list-style-type: none"> • Check all screws, connectors & parts are tightly fitted • Check all moving parts move freely • Check internal heating element connections are tight • Replace heating element if covered with lime scale • If plug, cable or socket are damaged, replace • Take care not to damage heater and sensors
Function Checks	<ul style="list-style-type: none"> • Inspect for defects e.g. cracks, severe pits in the chamber, kinks on tubes. All wires. • Clean and inspect water level sensor, door gasket • Test & inspect all valves e.g. fill, vent, air, repair or replace as needed • Check for leaks during operation • Check for operation of pressure gauges, temp. gauges, LED displays, LCD, and indicator lamps • Leakage current meets the safety standards • Earth lead is tightly secured to equipment

Troubleshooting

	Fault	Possible causes	Solution
1	Equipment NOT heating	No power at mains socket Electric cable fault Damaged heating element	Check power switch is on. Replace fuse. Check mains power socket is working. Rewiring if power not present. Try cable on another piece of equipment. Replace if broken
2	Pressure rises above the marked level	Blocked valve	Clean the pressure regulating valve, safety valve. Pressure vessel may be over filled Retest autoclave under pressure with water only.
3	Steam is constantly escaping	Poor seal	Clean leaky valve and hole, replace if defective. Clean leaking seal or gasket, replace if broken.
4	Electric shocks	Wiring fault	Rewiring/repair by electrician immediately

19) Hot Air Oven

Function

Hot Air Ovens have several functions depending on the set temperature.

In the Laboratory, they can incubate specimens as well as drying slides if set at low temperatures.

Sterilization of water sensitive items like powders, Vaseline, glassware as well as metallic instruments is possible at higher temperatures

How it works

Ovens are basically insulated enclosures in a metal box. They are electrically heated with an air circulating fan. A Thermostat up to 300 centigrade is fitted. A ventilation hole, a timer and a thermometer are also fitted. For the Lab, the thermometer controls the required temperature, with the ventilation hole open to allow steam escape.

During sterilization, the temperature is set high for a longer time to allow complete destruction of bacteria.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none">• Disconnect from power supply prior to cleaning• Clean the exterior case and interior chamber with a moist cleaning cloth• Dry the equipment thoroughly• Remove shelves and immerse in water and clean completely, and then dry shelves thoroughly Never use benzene or paint thinner for cleaning
Visual checks	<ul style="list-style-type: none">• Use a properly grounded electrical outlet• Ensure that outside of unit and power cord are clean and not damaged• Ensure that internal chamber is clean

	<ul style="list-style-type: none"> • Ensure that the Heater element, air inlets, vent outlets, control enclosures, door latches and hinges are clean
Function Checks	<ul style="list-style-type: none"> • Check that Temperature control sensor is not damaged • Check that circuit breakers & switches are operational • Check predetermined setting and temperatures • Significantly higher temperature could indicate temperature controller is failing

Troubleshooting

	Fault	Possible causes	Solution
1	Erratic temperature	Defective circuit board Dirty sensor	Repair or replace new one Clean the sensor
2	Can NOT heat	Defective circuit board Burnout heater Defective triac Air leakage	Replace new one Replace heater Replace triac Lock door completely or replace door gasket
3	Can NOT power up	Disconnected linkage Mal-functional switch Burnout heater Blew fuse	Connect linkage completely Replace switch Replace heater Replace fuse
4	Failure set temperature	Defective circuit board	Repair or replace new one
5	Unacceptable uniformity	Defective fan	Repair or replace fan or adjust the air vents
6	Temperature fluctuated	Dirty sensor, Burnout electronic parts Defective circuit board	Clean sensor Repair or replace parts Repair or replace circuit board

20) Disinfector/ Boiler

Function

Disinfection kills and reduces some microorganisms but not the resistant bacterial spores. A boiler is used to disinfect medical objects that come in contact with intact skin but not mucous membranes (none critical items). Examples of these items are; suction tubes, dishes and some instruments. Boilers can be electrically operated or non-electric (external heat source).

How it works

It has different parts that help it perform its work effectively. They include; chamber, cover with handle, equipment tray with handles, heating elements, electrical cable (electrical boiler), drain tap, thermal switch (for safety), thermostat (temperature control) and rubber stands.

Items to be disinfected first undergo a process of washing, cleaning and rinsing.

Clean items are then immersed into equipment tray, dipped in and properly covered under water and closed. Heat is introduced and water is left to boil at 100 degrees centigrade for twenty (20) minutes before it is switched off. The equipment tray is removed and water is drained out to allow for a fresh procedure to take place.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none">• After every use, wash and clean with vim.• Scrub scale lime deposits and wash with vim.• Change water every after use.• Clean and dry boiler when not in use.
Visual checks	<ul style="list-style-type: none">• Check for scaling on element and chamber.• Check for loose nuts on cover and drain tap. Check firmness of electrical connections of top plug and socket.

Function Checks	<ul style="list-style-type: none"> • Use trouble shooting guide (step by step approach). • Ensure thermostat operates and element heats. • Ensure switches are functional. • Check functionality before next use.
-----------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Troubleshooting

	Fault	Possible causes	Solution
1	Equipment NOT heating	No power at mains socket. Fault on electric cable. Heating element blown. Thermal safety switch tripped.	Replace fuse. Replace power cable. Replace heating element. Reset circuit breaker, switch on mains. Reset safety switch.
2	Water takes long to boil	One element faulty (for more than one element). Scale deposits on the element. Thermostat wrongly set.	Check one of the elements, replace. Remove scale deposits. Adjust thermostat setting.
3	Heating elements keep blowing frequently.	Much scale deposits. Little water in the chamber. Faulty safety switch. Faulty thermostat	De-scale heating element. Always put enough water/user train. Replace safety switch. Replace thermostat.
4	Much steam escapes.	Lid/cover does not close well. Drain valve does not close well. Many items in the boiler chamber	Tighten loose nuts and replace missing ones. Replace drain valve. Reduce on the items/train user.

21) ECG (Electrocardiogram) Machine

ECG machines are used to monitor the electrical activity of the heart and display it on a small screen or record it on a piece of paper. The recordings are used to diagnose the condition of the heart muscle and its nerve system.

The electrical activity is picked up by means of electrodes placed on the skin. The signal is amplified, processed if necessary and then ECG tracings displayed and printed. Some ECG machines also provide preliminary interpretation of ECG recordings. There are 12 different types of recording displayed depending upon the points from where the recordings are taken

Care must be taken to make the electrode sites clean of dirt before applying electrode jelly. Most problems occur with the patient cables or electrodes.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none">• Clean the Printing head• Clean body of equipment with a damp cloth and dry off
Visual Checks	<ul style="list-style-type: none">• Check that battery charge indicator, power indicator and patient cable connector indicators are working• Check that cables are not bent, knotted or damaged

	<ul style="list-style-type: none"> • Replace any damaged electrical plugs, sockets or cables • Check all knobs, switches and indicators are tightly fitted
Function Checks	<ul style="list-style-type: none"> • Calibrate machine before use using 1 mV pulse • Check that the baseline of the ECG recording is steady • Check the calibration of recordings with ECG a simulator • Check battery power can operate the equipment

Troubleshooting

	Fault	Possible causes	Solution
1	ECG traces have artifacts or base line drift	Improper grounding Patient is tensed up Deep breathing	Try with battery power only. If the recording improves then problem is with grounding. Power the machine from another outlet with proper electrical ground. Relax the patient
2	ECG traces have artifacts in one or more traces, but NOT in all traces	Improper electrode connection with patient or problem with the ECG cable	Check the patient cable. Replace cable if found faulty. Check the electrodes expiration. Check patient skin preparation. Check limb electrodes and chest electrodes for damage
3	Paper feed NOT advancing	Incorrect paper loading	Use instructions to reload paper
4	Printing NOT clear or NOT uniform	Printing head problem	Adjust the printing head temperature or position. Clean the printing head with head cleaner. If it does not improve, replace printing head. Check the paper roller and replace if not smooth
5	The machine shuts down after a few minutes while on battery power	Problem with battery or charging circuit	Recharge the unit overnight. If there is no improvement then replace the battery. If there is still no improvement, refer to a Technician

22) Pulse Oximeter

Function

A pulse oximeter is a device that non-invasively monitors the oxygen saturation of a patient's blood. It measures the amount of oxygen in a patient's arterial blood during operations and diagnosis. This level of oxygen or oxygen saturation is often referred to as SpO_2 , measured in %, and this is displayed on the pulse oximeter. A pulse oximeter also displays pulse rate.

How it works

The coloured substance in blood, haemoglobin, is a carrier of oxygen and the absorption of light by haemoglobin varies with the amount of oxygenation. Two different kinds of light (one visible, one invisible) are directed through the skin from one side of a probe, and the amount transmitted is measured on the other side. The machine converts the ratio of transmission of the two kinds of light into a % oxygenation. Pulse oximeter probes can be mounted on the finger or ear lobe.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none">• Clean the probe with alcohol wipe after each use• Unplug, clean outside with a damp cloth and dry off
Visual Checks	<ul style="list-style-type: none">• Check that all parts are present and connected• Check that cables are not twisted and remove from service if any damage is visible• Tighten any loose parts• If plug, cable or socket are damaged, replace
Function Checks	<ul style="list-style-type: none">• Check operation of all lights, indicators and visual displays• Check probe disconnection alarm.

Troubleshooting

	Fault	Possible causes	Solution
1	Equipment is NOT running	No power from mains socket Battery (if present) is discharged Electrical cable fault	Check power switch is on. Replace fuse. Check mains power is present at socket using equipment known to be working. Recharge or replace battery Try cable on another piece of equipment.
2	SpO ₂ or pulse rate NOT displayed or unstable	Probe is not mounted correctly Probe not able to read through dirt, nail polish, etc. Patient movement Patient's SpO ₂ value is too low to be measured	Connect probe and cable properly Remove grease, dirt, nail polish and clean probe Request patient to remain still Further clinical examination of patient. Recite probe if necessary
3	Probe off displayed on screen	Probe is not connected properly The connection between the probe and oximeter is loose	Connect the sensor Exchange cable
4	"Error" displayed on screen	Faulty probe or control circuit	Refer to user manual and troubleshoot by technician
5	Continuous alarm sounds	Alarm limits set too low or high Power disconnected	Set appropriate alarm limits Connect power cable

23) Slit Lamp

Function

A slit Lamp is used for observing the patient's eye.

The eyelid and anterior segment (cornea, iris, ciliary body, anterior chamber, posterior chamber, crystalline lens, lens, zinn's zonule, anterior vitreous, anterior sclera) is observed in the low magnification.


Posterior segment (posterior sclera, retina, choroid, posterior sclera) is observed in high magnification with Ruby lens.

How it works

By changing the slit width control knob or turning the aperture and slit length control wheel, ideal slit width, slit length and spot size can be achieved.

By turning the illumination unit or operating slit rotation control ring, the desired position of the slit image can be obtained.

Before use, diopter compensation and interpupillary distance adjustments are carried out. A focusing test rod is used to establish the proper microscope setting for each use. The test rod is inserted in the hole in chin-rest with the flat surface at the top of test rod facing the microscope.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none"> • Remove dust from the mirror and eye piece lens using camel hair brush. Finger marks can be removed if necessary with minimum isopropyl alcohol and a cotton swab. Dry with a facial tissue. • Clean exterior surfaces especially the joystick glide plate by wiping with a soft, dry cloth. Do not use commercial or household cleaners.
Visual checks	<ul style="list-style-type: none"> • If the glass envelope is blackening, replace the Halogen lamp. .
Function Checks	<ul style="list-style-type: none"> • Slit width control adjustment – loosen or tighten the small screw at the centre of the control to change the amount of force required to turn the knob. Should be tight enough to ensure that the slit does not close spontaneously. • Check the applanation tonometer at 0 g, 2 g and 6 g positions. • See the cone prism through left eyepiece and check whether cone prism is located at the centre. If not, adjust the cone position of vertical and horizontal direction. • Slit lamp blades should close without overlapping and open with parallel edges.

Troubleshooting

	Fault	Possible causes	Solution
1	No illumination.	Power cable, switch, fuse, bulb socket or bulb is damaged.	Check and replace damaged part.
2	Slit light is too dim	Bulb is not correctly inserted or Voltage selector setting is incorrect.	Insert bulb correctly, set Voltage selector to the correct position.

24) Centrifuge

Function

Centrifuges are used to fasten the sedimentation of substances in liquids. In medical laboratories specimen may be blood, urine etc. In the centrifuge, test tubes containing specimen are placed into buckets which are then rotated at high speeds or RPM (rotations per minute) for a pre-determined time. The substances in the test tubes are deposited in the order of weight, the heaviest element being the first to settle.

How it works

The brackets holding the buckets are coupled to the shaft of a motor which then rotates the shaft. The motor is normally a universal one with a variable speed mechanism. A timer is also incorporated in the circuit which cuts out power to the motor after a pre-set time by the user elapses. The two parameters being variable make it possible for various time/speed combinations to be applied to a varied range of specimen.

A safety feature is incorporated by the use of door switches. The switch completes the circuit only when the lid is closed to protect fingers from injuries that can be caused by the fast spinning head of the centrifuge.

Some larger centrifuges are equipped with an electrical or mechanical braking system for halting rotation'

Magnetic door locks in conjunction with centrifugal switches ensure that the lid cannot be opened as long as the head is rotating.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none"> • Clean the interior of the bowl with disinfectant daily.
Three monthly checks	<ul style="list-style-type: none"> • Check mountings and replace if necessary • Check for loose bolts and tighten • Check brushes and replace if necessary • Check for corrosion and repaint if necessary • Clean commutator with smooth sand paper

Troubleshooting

	Fault	Possible causes	Solution
1	Powers up but does NOT spin	<ul style="list-style-type: none"> • Worn out brushes • Dirty commutator 	Replace batteries Clean with smooth sand paper and contact spray
		<ul style="list-style-type: none"> • Check power supply and fuses 	Rectify accordingly
2	Noisy vibration	<ul style="list-style-type: none"> • Imbalance by missing buckets • Worn out bearings 	Ensure all buckets are in place Check and replace
3	Speed variation NOT responding	<ul style="list-style-type: none"> • Check the variable Resistor contacts 	Clean or replace

25) Microscope

Function

Microscopes are indispensable in medical laboratories. They enable the investigation of specimen (blood, urine etc.) for the presence of parasites or other abnormalities.

The microscope is able to achieve this by producing magnified images of specimen examined through it by means of multiple lens arrangement.

How it works

A microscope consists of two lenses (the objective and the eye piece) and a light source with a condenser. Rays from a light source are directed into a condenser which then brings them to a common focus on the specimen. The light rays then pass through the objective and produce a primary image in the focal plane of the eye piece. The eye piece magnifies the primary image and brings it into focus with the retina of the eye.

The objective and eye piece are systems of lenses where the distance from each other is fixed by the mechanical length of the tube. To adjust to form a clear image, either the specimen stage or tube as a whole has to be moved by means of coarse and fine adjustment drives.

The light source of a microscope can be either natural or artificial light.



Maintenance Checklist

Storage	<ul style="list-style-type: none"> • Should be kept in a dry room at temperature between 5°C to 35°C After use cover it with dust covers
Care and Handling	<ul style="list-style-type: none"> • Dust can be brushed off with a small paint bush • Clean dirty lenses with lens tissue or soft clean cloth. Breathe on the lens before wiping it. • If the pollution of the lenses is heavy use 2 to 3 drops of xylene on tissue paper • Check all cleaned lenses with a magnifier (inverted eye piece) • Clean oil immersion objective immediately after use with absorbent paper • Clean the body and housing with a soft cloth • Apply petroleum jelly on the rack and pinion when required

Troubleshooting

	Fault	Possible causes	Solution
1	Light source bulb NOT lighting	<ul style="list-style-type: none"> • Check the bulb • Intensity knob stuck on minimum 	Replace if blown Check the knob and rod assembly and clean
		<ul style="list-style-type: none"> • Check the light circuit fuses 	Replace if blown
		<ul style="list-style-type: none"> • Check the integrity of the transformer 	Replace if defective
2	Descending stage	Loose or dirty pinion and rack	Clean, lubricate and tighten the rack
		Check bearings between the stage and limb	Clean and lubricate

26) Laboratory Incubator

Function

Incubator is a device used to grow and maintain microbiological cultures or cell cultures. The incubator maintains optimal temperature, humidity and other conditions such as the carbon dioxide (CO₂) and oxygen content of the atmosphere inside.

The scope of usage of Incubators may vary but the basic requirement of incubator is to provide a controlled temperature environment.

How it works

Laboratory incubators consist of a chamber, heat source, a resettable thermostat and a thermometer. Heat generated is controlled by the thermostat to a desired temperature that is set by the user. The thermometer displays the temperature in the chamber.

Maintenance Checklist

Care and Handling	<ul style="list-style-type: none">• Use a smooth cloth soaked in water to clean the surface of the equipment.• Disinfect the machine with an appropriate disinfection solution e.g. Hypochlorite (JIK)
Visual checks	<ul style="list-style-type: none">• Make sure the equipment shows no error messages on the controls.• Check whether the door locks properly.
Function Checks	<ul style="list-style-type: none">• Check whether the whole system functions before use (e.g. temperature control and display, power supply).

Troubleshooting

	Fault	Possible causes	Solution
1	No display on control panel	No power	The power cord is loose in the socket outlet. Check the outlet for power. Check the circuit breaker on the power panel. In case of tripping, reset it.

			Check that the main power switch on the incubator is turned on. Check for any blown fuses
2	High temperature alarm flashes	Faulty thermostat	Check if temperature probe is well connected Check whether the thermostat is able to switch on and off
3	Oxidation forming on the interior surfaces	Use of high concentrated detergent.	Clean the interior surface with soft water
4	Temperature display is drifting	Location effect	Check the incubator location.
5	Inside temperature is NOT equal to the display temperature.	Needs calibration Location effect	Verify the calibration procedure. Check the location of the incubator.
6	Temperature NOT controlling at set point	Fault with set point Chambers circulated with cold air	Check the set point Close the doors for a minimum of 30–40 minutes
7	Display on but NO temperature rise	Fault on the Thermostat Fault on the Element	Check thermostat whether it's okay if not replace Check element, if blown replace
8	Temperatures take long to reach the set temperature	Faulty door system	Check the gasket, if damaged replace.

27) Colorimeter

Function

A colorimeter is used to measure the concentration of a substance in the patient's sample by comparing the amount of light it absorbs with that absorbed by a standard sample that contains a known amount of the substance being tested.

A coloured solution of the substance being measured or a coloured derivative of it is produced. Coloured solutions absorb light at given wavelengths in the visible spectrum.


How it works

The essential parts of a colorimeter are a light source, an adjustable aperture, colored filters, a cuvette holder, a detector and a data display.

Colored filters are used to select the required wavelength. The usual wavelength range is from 400 to 700 nanometers (nm).

A cuvette which obtains light path of 10 mm size is recommended. The data display indicates the results as transmittance (a linear scale from 0–100%) or as absorbance (a logarithmic scale from zero to infinity).

To prolong the life of the lamp, switch off the device after use and disconnect from the mains socket at the end of each day.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none"> • Clean the equipment with a soft cloth soaked in a mild disinfection reagent • Clean the cuvette holder
Visual checks	<ul style="list-style-type: none"> • Check the functionality of the bulbs. • Check that there is no mechanical damage • Check that the cuvette is free of scratches and finger marks
Function Checks	<ul style="list-style-type: none"> • Switch ON/OFF and change the filters smoothly • Check that the correct type of cuvette is in use • Check that the display is functional when switched ON

Troubleshooting

	Fault	Possible causes	Solution
1	The machine does NOT get power	Faulty power out let socket Power cord loosely connected Blown fuses	Check the power socket, if faulty, replace. Fix the power cord. Replace fuses of the same rating.
2	The machine powers but no light	Blown bulb	Replace the bulb of the same rating.
3	The machine is NOT zeroing	Machine is due for calibration Dust accumulation on the light filter	Carry out calibration with a standard medium Clean the filters.

28) Haematology Analyzer

Function

A Haematology analyzer is an automatic multi parameter blood cell counter for in vitro diagnostic use in clinical laboratories. It performs speedy and accurate analysis of parameters in blood and detects the abnormal samples. The instrument displays abnormal analysis data with abnormal marks attached on the LCD screen. Displayed analysis data allows detecting those samples which are outside the tolerance and need further analysis and consideration. Whole blood and pre-diluted blood can be analyzed.

How it works

The white blood cell (WBC) count is measured by the WBC detector block using DC detection method. The red blood cell (RBC) count and platelets are taken by the RBC detector block, also using the DC detection method. The HGB detector block measures the hemoglobin concentration using the non-cyanide hemoglobin method.



Maintenance Checklist

Caring and Handling	<ul style="list-style-type: none">• Clean the transducer (TD) chamber and diluted sample line daily(Execute shut down).• Clean the sample rotary valve (SRV) tray weekly.• Clean the waste chamber (rinse sequence) monthly or after every 2500 samples.• Clean the transducer (rinse sequence) monthly or after every 2500 samples.• Clean the sample rotary valve (SRV) every 3 months or after 7500 samples. <i>SRV</i> life time is 5 years• Clean the rinse cup as needed.
---------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

	<ul style="list-style-type: none"> • Clean the WBC/RBC transducer aperture as needed.
Visual checks	<ul style="list-style-type: none"> • Check trap chamber level daily and discard if necessary. • Auto rinse as needed. • Replace the waste tank as needed.
Function Checks	<ul style="list-style-type: none"> • Perform a reagent replacement sequence • Auto rinse with background check sequence • Check settings sub-menu for system set up • Perform the waste chamber cleaning sequence • Perform the transducer cleaning sequence • Perform transducer fluid draining sequence (for clog removal) • Check status display <ol style="list-style-type: none"> 1. Hemoglobin (HGB) convert (real time) 2. Pressure and Vacuum (real time) 3. Unit operation counter • Check paper feed (optional built-in printer)

Troubleshooting

	Fault	Possible causes	Solution
1	Alarm sounds and an error message is displayed.	Pressure or vacuum, is outside range. Air bubbles in the system. Clogging in transducer aperture(s)	Press [HELP] key on keyboard to stop the alarm and change over to the HELP screen that shows what actions to take against the error.

29) CD4 Counter, Flow Cytometer Type

Function

CD4 counter uses flow Cytometry technology that simultaneously measures/counts and analyses multiple characteristics of single cells or particles such as CD4 (One of the lymphocyte). The CD4 count data provides important information for staging and monitoring patients infected with HIV/AIDS. The software provides an automated analysis without operator's intervention. The normal CD4 count is 1,000 cells /mm³ and this count is lower for HIV infected patients.

How it works

The tissue sample is broken up into single cells and held in a test tube, which is placed into the flow cytometer. The liquid containing the cells is drawn up from the test tube and pumped to the flow chamber (flow cell). Cells flow through the flow chamber one at a time very quickly and are presented to one or more light sources (Lasers).

A small laser beam of very bright light hits the cells as they pass through the flow chamber. The pattern of light scattering is dependent on cell size and shape giving relative measurements of these cellular morphological characteristics as cells flow through the beam. Forward scatter measures light scattered in the direction of the laser path and measures the size of the cell. Side scatter measures scattered light at 90 degrees to the laser path and measures the cell granularity.

The light detector processes the light signals and sends the information to the computer. Each type of cell in the immune system has a unique combination of forward and side scatter measurements, allowing count of number of each type of cell. Filters are used to direct the light emitted by the fluorochromes to the colour detectors. The fluorochromes attached to the cells absorb light and then emit a specific colour of light depending on

the type of fluorochrome. The colour detectors (PMT's) collect the different colours of light emitted by the fluorochromes.

The electronics in the cytometer amplify and process the resulting data. They convert analogue data to digital data which is stored in the computer.

This data from the light detector and the colour detectors is analysed from the computer and plotted on histograms.



Maintenance Checklist

Care and Handling	<ul style="list-style-type: none"> • Unplug, clean external surfaces with a damp cloth with cleaning solution containing alcohol for example 70% ethanol, Deionized water and wipe surfaces with a clean, dry cloth
Visual checks	<ul style="list-style-type: none"> • Check that all parts are there and correctly fitted • Check that the mains cable has no damage and if damaged replace • Check that the air filter is free of dust and if clogged wash clean with water, dry and reinstall or replace • Check the level of sheath fluid in the reservoir • Check the level of waste in waste reservoir , dispose of if necessary
Function Checks	<ul style="list-style-type: none"> • Check photo multiplier tube (PMT) voltage adjustment • Check time – delay calibration • Check fluorescence compensation adjustment • Carry out sensitivity testing

Troubleshooting

	Fault	Possible causes	Solution
1	Flow cell NOT filling	No sheath pressure	Check for leaks or cracks in the sheath reservoir. Replace it, if needed. Tighten the sheath cap. Verify the vent valve is in forward
2	No events displayed and status = READY	Communication failure between computer and instrument.	Turn off the computer and instrument; Turn on the instrument, then the computer.
3	No events displayed and status = READY	No sample in the sampler tube	Add sample or install new sample tube
4	Scatter parameters distorted	Air bubble in flow cell Incorrect instrument settings	Prime the fluidics Press PRIME Rerun BD FACSComp
5	Time –delay calibration fails	Empty sheath reservoir. Air in sheath filter Air bubble in the flow cell	Fill the sheath reservoir Vent air from the sheath filter Prime the fluidics
6	Sample tube NOT pressurizing	Fluid mode not in RUN Cracked sample tube Worn Loader seal	Press RUN Replace the sample tube Replace the seal

30) Dental Unit and Dental Chair

Function

Dental unit is intended for diagnostic and therapeutic treatment of dental patients by dentists or health care professionals. Generally, a dental chair refers to the chair for patients while the dental unit consists of the chair, treatment hand pieces, a doctor's table, a light and a cuspidor unit.

How it works

1. Turn the Master ON/OFF toggle to the "ON" position.
2. The hand pieces are automatically activated when removed from their hanger.
3. Depress the Foot Control to operate the activated hand piece. Pressure is shown on the gauge for the hand piece being used.
4. To adjust the pressure to each hand piece, turn the adjustment screw clockwise to decrease pressure or counter-clockwise to increase pressure.

CAUTION: When adjusting the hand piece pressure, do not over-tighten the screws.

5. To activate water spray, turn the water ON/OFF Toggle to the "ON" position. Flow adjustment to the "wet" hand pieces is controlled by the controls labelled "Flow Control".
6. To release the air of the unit arm, activate the momentary toggle valve and adjust the height of the flex arm to desired position. When desired position is achieved, release the momentary toggle switch to lock the height of the flex arm.

CAUTION: Do not attempt to adjust the height of the flex are without releasing the air-brake valve. Failure to release the air-brake valve may cause damage.



Maintenance Checklist

Cleaning [Daily]	<ul style="list-style-type: none"> ● Purge the unit with air ● Remove hand pieces from the tubing. ● Empty the water bottle, then reinstall it (if water bottle is applicable, when using city water, turn off water from the source) ● Hold the hand piece tubing and syringe over a pail. Turn the unit on, wait a moment and then operate the flush toggle, syringe and foot control unit water is purged from the system ● Turn unit off ● Disinfecting the bottle; Fill the bottle with 100ml disinfectant solution (9 parts tap water & 1 part 5.25% Sodium hypochlorite/ household bleach), shake vigorously and let settle for 10 minutes. Shake again and then rinse twice with water.
Cleaning [Weekly]	<ul style="list-style-type: none"> ● Purge the unit with air ● Flush the system with disinfectant solution <ol style="list-style-type: none"> a. Turn unit off. Empty the water bottle, replacing the water with cleaning solution b. Remove hand pieces from tubing and hold the hand piece tubing and syringe over a pail. c. Turn unit on, wait a few moments and then operate the flush toggle, syringe and foot control unit a continuous stream of solution is running through the system ● Allow the disinfectant to remain in the unit for at least 10 to 20 minutes and then flush the system again unit all the solution is used up. ● Purge the unit with air; <ol style="list-style-type: none"> a. Hold the hand piece tubing and syringe over a pail. Turn the unit on, with a few moments and then operate flush toggle, syringe and foot control unit all solution is purged from the system. b. Turn unit off. ● Fill with clean water ● Turned unit off, remove the empty disinfectant bottle. Replace with clean bottle and water.

	<ul style="list-style-type: none"> Disinfectant solution; Use 100ml of disinfectant solution.
Visual checks	<ul style="list-style-type: none"> When the compressor is OK, we expect lighting system to be OK.
Function Checks	<p>When the power is on you check</p> <ul style="list-style-type: none"> Hand pieces, Air compressor, Foot switch Check knobs are operational

Troubleshooting

	Fault	Possible causes
1	Hand piece lacks Power	Check regulator adjustment (80psi) Check hand piece pressure adjustment on control block Pinched supply Bad hand piece gasket at connection with tubing Defective hand piece
2	Water coolant does NOT shut off when foot switch is released	Adjust air pressure to 80psi – Water pressure to 40psi Foot control is not exhausting Defective water relay in valve
3	More than one hand piece is operating	Hand piece is not completely in the hanger Improve adjustment of pilot valve in the hanger. Kinked or pinched signal line from the pilot valve.
4	Insufficient water coolant	Adjust coolant flow valve Water filter may be plugged Plugged hand piece Kinked or pinched tubing Improper adjustment of water relay
5	Water coolant is running from hand piece while in holder	Water pressure is too high Air pressure is too low Hand piece holder out of adjustment
6	Water coolant is running continuously	Purge switch is on Water pressure is too high Hand piece holder out of adjustment Improper adjustment of water relay

31) Solar Photovoltaic (PV) System

Function

Solar PV systems generate electricity to operate recommended appliances (e.g. lamps for lighting, ultrasound scanner, oxygen concentrator, suction machine, computer, etc.)

How it works

A solar PV system uses sunlight to generate electricity. A solar PV system consists of the following main components:

Solar panel/module:

- It receives sunlight and converts it into electricity to charge the battery.

Charge controller (Regulator):

- It controls battery charging and discharging.
- It prevents overcharge (too much current into the battery) of the battery by automatically disconnecting the solar panels when the battery is fully charged.
- It disconnects loads to prevent battery discharge (draining of the battery) beyond safe levels.

Battery:

- It stores electricity generated by the solar panels for use when required.

Inverter:

- It converts the direct current (DC) from the battery into alternating current (AC) and supplies it to operate AC appliances.

AC and DC loads (Appliances):

- These are the AC and DC loads that use electricity.

Other Accessories (Distribution board, Circuit breakers, Cables):

- These are used to interconnect, control and protect different components against damage.

Maintenance Checklist

Care and Handling	<ul style="list-style-type: none"> • Clean the surface of the panels with water and soft cloth to remove any dirt. • Clean the battery terminals of any dust, sulphate deposits and apply a thin layer of Vaseline to avoid corrosion. • Ensure that all the lights are switched off every morning.
Visual checks	<ul style="list-style-type: none"> • Check that all LED lights are functional • Check that the solar modules are not cracked or shaded. •
Function Checks	<ul style="list-style-type: none"> • Check that the LED lights (system operation information light and state of charge on battery light) in the morning and evening are green.


Troubleshooting

	Fault	Possible causes	Solution
1	System operation status LED light in Red.	<ul style="list-style-type: none"> • Fuse is blown. • System not charging or regulator has over heated. 	<ul style="list-style-type: none"> • Replace the fuse. • Check the regulator functionality and ensure solar panels are not shaded or clogged by dirt. •
2	Battery status light in Red	<ul style="list-style-type: none"> • Low battery charge. • Charge regulator short circuiting 	<ul style="list-style-type: none"> • Disconnect loads and allow the battery enough time to charge and if it fails, replace the battery. • Replace the charge regulator.
3	Both LED lights are green but no power output.	<ul style="list-style-type: none"> • The inverter could be off. • The circuit breaker could have tripped. 	<ul style="list-style-type: none"> • Switch on the inverter. • Put back the breaker to the on position. • Trace for overload and disconnect.

Annexes

- Annex 1: Standard Floor Plan of a Regional Workshop
- Annex 2: Recommended RW Tools, Test Equipment and Furniture
- Annex 3: Equipment Condition A–F
- Annex 4: Medical Equipment Inventory Form
- Annex 5: Annual Work Plan Form
- Annex 6(a): 5S Activity's Instruction for RWs
- Annex 6(b): CQI Action Plan Form
- Annex 7: Quarterly Report Form
- Annex 8: Job Card Form
- Annex 9: Annual Assessment Sheet for RWs
- Annex 10: Maintenance Requisition Form
- Annex 11: Maintenance Service Sticker Form

Annex 1: Standard Floor Plan of a Regional Workshop



Annex 2: Recommended RW Tools, Test Equipment and Furniture

No.	Tools and Test Equipment	Qt'y
1	Tool Kit, Mechanical ^(*)	2
2	Tool Kit, Electrical ^(*)	2
3	Tool Kit, Electronics ^(*)	2
4	Tool Kit, Biomedical ^(*)	2
5	Tool Kit, Plumber ^(*)	2
6	Tool Kit, Refrigeration ^(*)	2
7	Tool Kit, Carpentry	1
8	Tool Kit, Civil	1
9	Biosafety cabinet testing and calibration toolkit	1
10	Welding machine, Heavy Duty	1
11	Welding machine, portable	1
12	Gas welding machine	1
13	Workshop stool	4
14	Phototherapy Light Intensity (Spectral Irradiance) Meter	1
15	Soldering Kit, 2 in 1 Rework soldering station with hot air and heat gun.	1
16	De-soldering kit, standard	1
17	Spanner Set (Metric)	1
18	Spanner Set (Imperial)	1
19	3-piece extendable ladder/foldable ladder	1
20	Air Compressor	1
21	Spray gun	1
22	AC/DC regulated power supply, DC to DC convertor	1
23	Battery Charger	1
24	Hydrometer	1
25	Battery tester	1
26	Threading toolkit.	1
27	Ventilator/Anaesthesia analyzer	1
28	Workshop protective wear (helmets, eyewear, gloves, shoes, aspirators overcoat/overall/trouser and shirt)	10
29	Automated electrical safety analyzer	1
30	Ultrasound electrical safety transducer leakage current tester	1
31	Portable oxygen analyzer	1
32	NIBP, SpO ₂ , ECG simulator	1

No.	Tools and Test Equipment	Qt'y
33	Infusion pump analyzer	1
34	ESU (Electro–Surgical Unit) Analyzer	1
35	Defibrillator Analyzer	1
36	mA meter	1
37	Digital manometer	1
38	2 Drawers work bench with vice No. 5	1
39	Engraving machine, electric, heavy duty	1
40	13mm chuck, reversible, variable speed, Drilling machine, with hammer function, electrical, hand held	1
41	13mm chuck, reversible, 2 speed, Portal drilling machine, cordless with spare battery and charger	1
42	Vice No. 3	1
43	Portable Generator set, 5KVA minimum	1
44	Portable Grinder	1
45	Dual DC regulated power supply	1
46	3 Phase sequence tester	1
47	Mega tester	1
Office Furniture		
48	Office Chair	2
49	Chair	10
50	Desk with lockable drawers	3
51	3 in 1 printer/scanner/photocopier	1
52	Filing cabinet	3
53	Desktop Computer	2
54	Reliable Internet Connectivity	1
55	Phoneline connection to all wards	1
56	6 Persons conference table	1
57	Shelves	3
58	White board 1600cm x 110cm	1

Note: (*) The details of the Kits are described in the following pages.

No. 1/Tool Kit, Mechanical

No.	Description (Tool Kit, Mechanical)	Qt'y
1	Cutter, side, large	1
2	G-clamp	1
3	Pliers, long nose, 170mm	1
4	Screwdriver/tester, Voltage tester	1
5	Box, tool, steel, big, 4 tray, cantilever, 550 x 220 x 215mm (l x w x h)	1
6	Oil can, 125ml, metal with pump	1
7	Hydrometer/battery electrolyte tester	1
8	Brush, steel, hand, 290mm; 37mm; 4 rows	1
9	Chisel and punch set, 6 pc	1
10	File set flat/square/half round/round, coarse and fine, 200mm	1
11	Hacksaw, standard, 470mm	1
12	Hammer, ball pein, 500g/1lb; 350mm	1
13	Hammer, nylon, 280mm; diam:32mm	1
14	Key, hexagon, metric, set, 9pc, 1.5 – 10mm	1
15	Machine, drilling, hand, 2 speed, 1 – 12mm	1
16	Pliers, water pump, standard, 225mm	1
17	Punch, nail, 150 x 5 mm	1
18	Letter punch set, A-Z	1
19	Number punch set, 0-9	1
20	Screwdriver, flat, blade: 3.5mm; shaft:75mm	1
21	Screwdriver, flat, blade: 5.5mm; shaft 100mm	1
22	Screwdriver, flat, blade: 8mm; shaft: 250cm	1
23	Screwdriver, star, No. 1x, Diam: 4 mm; shaft 75mm	1
24	Screwdriver, star, No. 2x, Diam: 6mm; shaft 100mm	1
25	Screwdriver, star, No. 3x., Diam: 8mm; shaft: 150mm	1
26	Spanner, flat, metric, set, 15 pc, 6 – 32mm	1
27	Spanner, ring/flat, metric, set, 16pc, 5 – 25mm	1
28	Spanner, box/socket, metric set, 16 pcs, 6-33mm	1
29	Wrench, self gripping, Max: 50mm	1
30	Wrench, socket, set, large, 10 – 32mm; ½” square drive	1
31	Wrench, chain, diam.: 7 – 2200mm	1
32	Tap and die set, imperial	1
33	Tap and die set, metric	1
34	Wire brush	1
35	Reamer set	1
36	Steel rule, marked, metric: 300mm	1
37	Tape measure, Steel, 3m	1

No. 2/Tool Kit, Electrician

No.	Description (Tool Kit, Electrician)	Qt'y
1	Cutter, side, large, Max. size wire: 2.0mm; L = 210mm	1
2	Multi-meter, basic, digital, Like fluke 11	1
3	Bit, drill, concrete, hammer operated, 3.5/4.5/5.0/6.0 mm, with handle	1
4	Hammer, club, 42 x 42 x 110mm	1
5	Screwdriver, flat, insulated, blade 2.5mm; shaft 75mm	1
6	Screwdriver, flat, insulated, blade 4.0mm; shaft 100mm	1
7	Screwdriver, star, insulated, blade 2.5mm; shaft 75mm	1
8	Screwdriver, star, insulated, blade 4.0mm; shaft 100mm	1
9	Screwdriver/ mains tester, blade 4mm; shaft: 100mm, 100–500V	1
10	Stripper, Wire, Standard, 170mm	1
11	Knife, trimming, retractable blade	1
12	Pen, marking, waterproof, medium	1
13	Tape measure, Steel, 3m	1
14	Box, tool, steel, big, 4 tray, cantilever, 550 x 220 x 215mm (l x w x h)	1
15	Chisel, channeling	1
16	Chisel set complete with floorboard	1
17	Cutter, pipe, 3 – 30mm	1
18	File set flat/square/half round/round, coarse and fine, 200mm	1
19	Hacksaw, Junior	1
20	Hacksaw, standard, 470mm	1
21	Light, torch, rubber coated, 2 batteries LR20, 3V	1
22	Pliers, combination, 203mm / 8"	1
23	Pliers, long nose	1
24	Punch, centre, diam: 6mm; 1 = 115mm	1
25	Reamer, pipe	1
26	Screwdriver, flat, chubby, blade: 5.5mm; shaft: 40mm	1
27	Screwdriver, flat, long, blade: 5.5mm; shaft: 200mm	1
28	Screwdriver, flat, screw gripping, blade: 6mm; shaft: 200mm	1
29	Screwdriver, Phillips, chubby, No. 2x Diam: 5.0mm; shaft: 40mm	1
30	Screwdriver, Pozidriv, No. 1, 75mm	1
31	Screwdriver, Pozidriv, No. 2, 100mm	1
32	Spanner, adjustable, medium, 205mm / 8"; jaw:27m	1
33	Spectacles, safety	1
34	Wrench, pipe, small, 229mm / 9"	1
35	Malet Hammer	1
36	Claw Hammer	1
37	Wire brush	2
38	Grip pliers, adjustable, set	1
39	Hygrometer	1
40	Portable hand drill, cordless with rechargeable battery & charger	1
41	Allen keys set (metric & imperial)	1
42	Circlip plier set	1

No. 3/Tool Kit, Electronics

No.	Description (Tool Kit, Electronics)	Qt'y
1	Cutter, side, small size, L=125	1
2	De-soldering tool, standard	1
3	Dispenser, solder	1
4	Iron, soldering, medium, general, 30W / 220V; tip diam: 6mm	1
5	Multi-meter, digital, Like fluke 79 series, Resol: 3.5 digit	1
6	Pliers, long nose, 170mm	1
7	Screwdriver / mains tester, blade: 4mm; shaft: 100mm, 100 – 500V	1
8	Stripper, wire, standard, 170mm	1
9	Tools, trimming, set	1
10	Wick, de-solder	1
11	File, needle, set, 6pcs: 160mm	1
12	Knife, trimming, retractable blade	1
13	Tweezers, 130mm; tips: 1mm	1
14	Box, tool, brief-case type, 500 x 380w x 150h mm	1
15	File, half round, second cut, 150mm	1
16	Hacksaw, Junior	1
17	Hammer, ball pein, 100gr/1 oz	1
18	Key, hexagon, metric, set, 9pc, 1.5 – 10mm	1
19	Key, hexagon, set, imperial, 9pc, 5/64 – 3/8”	1
20	Magnifying glass, pocket type	1
21	Mirror, inspection	1
22	Pliers, combination; 180mm / 7”	1
23	Screwdriver, flat, blade: 3.5mm; shaft 75mm	1
24	Screwdriver, flat, blade: 5.5mm; shaft 100mm	1
25	Screwdriver, flat, blade: 6.5mm; shaft 100mm	1
26	Screwdriver, flat, long shaft, blade: 4.0mm; shaft 250mm	1
27	Screwdriver, flat, screw gripping, blade: 6mm; shaft 200mm	1
28	Screwdriver, Pozidriv. No. 0, 60mm	1
29	Screwdriver, Pozidriv. Long shaft, No. 1, 250mm	1
30	Screwdriver, Pozidriv. Long shaft, No. 2, 250mm	1
31	Screwdrivers, jeweller's, set, crosshead	1
32	Screwdrivers, jeweller's, set, flat blade	1
33	Shears, light duty, L 190 D 4.5	1
34	Spanner, adjustable, small, 155mm / 6”, jaw: 20mm	1
35	Spanner, ring/flat, metric, set, 8pc, 8 – 19mm	1
36	Brush, cleaning, 12mm, 260mm	1
37	Snipper set	1
38	Allen key set (metric & imperial)	1
39	Precision screw driver set	1

No. 4/Tool Kit, Medical Equipment

No.	Description (Tool Kit, Medical Equipment)	Qt'y
1	Cutter, side, small size, L=125	1
2	De-soldering tool, standard	1
3	Dispenser, solder	1
4	Iron, soldering, medium, general, 30W / 220V; tip diam: 6mm	1
5	Clamp meter	1
6	Pliers, long nose, 170mm	1
7	Screwdriver / mains tester, blade: 4mm; shaft: 100mm, 100 – 500V	1
8	Stripper, wire, standard, 170mm	1
9	Tools, trimming, set	1
10	File, needle, set, 6pcs: 160mm	1
11	Knife, trimming, retractable blade	1
12	Tweezers, 130mm; tips: 1mm	1
13	Box, tool, brief-case type, 500 x 380w x 150h mm	1
14	File, half round, second cut, 150mm	1
15	Hacksaw, Junior	1
16	Key, hexagon, metric, set, 9pc, 1.5 – 10mm	1
17	Magnifying glass, pocket type	1
18	Mirror, inspection	1
19	Pliers, combination; 180mm / 7"	1
20	Wrench, adjustable, medium, 350mm	1
21	Wrench, adjustable, small, 200mm	1
22	Wrench, socket, medium, 350mm	1
23	Wrench, socket, small, 200mm	1
24	Screwdriver, flat, blade: 3.5mm; shaft 75mm	1
25	Screwdriver, flat, blade: 5.5mm; shaft 100mm	1
26	Screwdriver, flat, blade: 6.5mm; shaft 100mm	1
27	Screwdriver, flat, long shaft, blade: 4.0mm; shaft 250mm	1
28	Screwdriver, flat, screw gripping, blade: 6mm; shaft 200mm	1
29	Screwdriver, Pozidriv. No. 0, 60mm	1
30	Screwdriver, Pozidriv. Long shaft, No. 1, 250mm	1
31	Screwdriver, Pozidriv. Long shaft, No. 2, 250mm	1
32	Screwdrivers, jeweller's, set, crosshead	1
33	Screwdrivers, jeweller's, set, flat blade	1
34	Spanner, adjustable, small, 155mm / 6", jaw: 20mm	1
35	Spanner, ring/flat, metric, set, 8pc, 8 – 19mm	1
36	Brush, cleaning, 12mm, 260mm	1
37	Snipper set	1
38	Allen key set, metric	1
39	Precision screw driver set	1

No. 5/Tool Kit, Plumber

No.	Description (Tool Kit Plumber)	Qt'y
1	Bit, drill, concrete, set, hammer operated, 3.5, 4.5, 5, 6mm; with handle	1
2	Level, Spirit, L = 50cm	1
3	Trowel, pointing, 150mm	1
4	Iron, brazing/soldering, large, 300W / 220V	1
5	Tape measure, steel, 3m	1
6	Box, tool, steel, big, 4 tray, cantilever, 550 x 220 x 215mm (l x w x h)	1
7	Chisel, cold, 16mm x 250mm	1
8	Chisel, floorboard	1
9	Cutter, pipe, 3 – 30mm	1
10	Hacksaw, junior	1
11	Hacksaw, standard, 470mm	1
12	Lamp, soldering, paraffin, 0.5 litre	1
13	Mallet, bossing, 60mm	1
14	Mallet, tinmans, 60mm	1
15	Pliers, combinations, 203mm / 8”	1
16	Pliers, water pump, standard, 225mm	1
17	Rule, s/s, 300mm 25mm wide	1
18	Shears, metal	1
19	Wrench, footprints, 220mm	1
20	Hammer, plumbers' 500gr	1
21	Hook, shave, triangular	1
22	Plumber's bob, with line, 100 gr	1
23	Portal Tripod stand	1
24	Pipe cutter	1
25	Wrench, adjustable, medium, 350mm	1
26	Wrench, adjustable, small, 200mm	1
27	Wrench, socket, medium, 350mm	1
28	Wrench, socket, small, 200mm	1
29	Bit, drill, auger, wood, set, 5pc, 6 – 25mm / 125mm	1
30	Brace, bit, ratchet, 4 jaw, 255mm	1

No. 6/Tool Kit, Refrigeration

No.	Description (Tool Kit, Refrigeration)	Qt'y
1	Leak detector, electronic	1
2	Pressure gauge (for R12 & R134a refrigerants)	1
3	Pipe bender set (1/2", 3/4", 1" ... 2")	1
4	Clamp meter	1
5	Flaring swaging tool	1
6	Ratchet spanner set	1
7	Portal brazing cylinder set	1
8	Blower (dust)	1
9	Digital thermometer	1
10	Portable brazing torch (butane)	1
11	Allen keys (metric & imperial)	1
12	Protective goggles	2
13	Portable gas welding kit	1
14	Silicon gun	1
15	Rivet gun	1
16	Tape measure, steel, 5 m	1
17	Electrical drill, with both hammer & reverse functions	1
18	Cordless portal hand drill, with rechargeable batteries	1
19	Digital clipper	1
20	Screw driver set (flat & screw)	1
21	Vacuum pump	1
22	Circlip pliers	1
23	Electrical clumping tool	1
24	Reamer set Refrigerant tester	1

Annex 3: Equipment Condition A-F

**CLASSIFICATION OF CONDITION
FOR MEDICAL EQUIPMENT
BY COLOR CODE**

 GOOD AND IN USE	NO PROBLEM
 GOOD BUT NOT IN USE	USER TRAINING NEEDED
 IN USE BUT NEEDS REPAIR	MAINTENANCE NEEDED
 IN USE BUT NEEDS REPLACEMENT (OLD OR OBSOLETE)	PLAN FOR DECOMMISSION -ING/DISPOSAL
 OUT OF ORDER BUT REPAIRABLE	MAINTENANCE NEEDED
 OUT OF ORDER AND SHOULD BE REPLACED (CAN NOT BE REPAIRED, NO SPARE, OLD OR OBSOLETE)	DISPOSE

Annex 5: Annual Work Plan Form

Output	Planned Output Target	Planned Activities to Deliver Outputs	Activity timing				GOU Quarterly allocation				(Name of IP) Quarterly allocation				TOTAL		
			Q1	Q2	Q3	Q4	Q1	Q2	Q3	Q4	Q1	Q2	Q3	Q4			
080803 (Maintenance of Medical and Solar Equipment)	Carry out PPM and PCM in 2RRHs, 4GHs and 9HCIVs	PPM and PCM in hospitals and HCIVs in ** region		8,000,000	12,000,000	12,000,000	8,000,000										40,000,000
080403 (Maintenance of Medical and Solar Equipment)	Procure assorted Medical Equipment spare parts	Procurement of assorted Medical Equipment spare parts		8,000,000	25,000,000		25,000,000										75,000,000
080803 (Maintenance of Medical and Solar Equipment)	Train 25 BME/Techs	Biomedical Engineering training in maintenance of ICU and theatre equipment.		12,000,000					12,000,000								24,000,000
080803 (Maintenance of Medical and Solar Equipment)	Collect and update equipment inventory database.	Collect and update of Medical equipment and solar systems inventory.		3,000,000			3,000,000										6,000,000
080801 (Technical support, monitoring and evaluation)	Conduct quarterly RW performance review meetings	Attend regional equipment Workshops' performance review meetings		2,400,000	2,400,000	2,400,000	2,400,000										14,400,000
080801 (Technical support, monitoring and evaluation)	At least 1 meeting organised and held	Hold regional medical equipment maintenance workshop management committee meetings					18,000,000										18,000,000
(Name of IP) operational budget support	Enhance technical capacity for Maintenance of lab equipment	Provide oversight support and mentorship											1,500,000				3,000,000
(Name of IP) operational budget support	NOMAD database system for equipment inventory management.	Upgrade NOMAD database to web based equipment inventory system and supervise its implementation.											3,000,000				3,000,000
TOTAL				33,400,000	39,400,000	35,400,000	35,400,000		13,200,000	17,700,000	4,200,000	2,700,000					181,400,000

Annex 6(a): 5S Activity's Instruction for RWs

1. What is 5S?

5S is a management tool, which originated from the Japanese manufacturing sector. It is used as a basic, fundamental, systematic approach for productivity, quality and safety improvement in all types of organizations.

5S is literally five abbreviations of Japanese terms with five initials of S. There are *Seiri*, *Seiton*, *Seiso*, *Seiketsu* and *Shitsuke*. In English, 5Ss were translated as *Sort*, *Set*, *Shine*, *Standardize* and *Sustain*. These are explained briefly below:

1. <u>Sort</u>	Remove unused stuff from your venue of work and reduce clutter	(Removal)
2. <u>Set</u>	Organize everything needed in proper order for easy operation	(Orderliness)
3. <u>Shine</u>	Maintain high standard of cleanness	(Cleanness)
4. <u>Standardize</u>	Set up the above three So as norms in every section of your place by use of Standard Operating Procedures and checklists	(Standardize)

The five steps of Sort–Set–Shine–Standardize–Sustain are a sequence of activities to improve the work environment to be as convenient and comfortable as possible and thereby also improve service contents with respect to preparedness, standardization and timeliness.

2. Practical hints for five steps

The basic actions and practical hint for implementing 5S are shown in the following page. The principal example and practical hints are shown in the following page for every step of 5S.

1S: SORTING / Elimination

Basic actions for sorting



1. Create a chaos free environment. 3. Use “Red Tags.”
2. Clean the walls and notice boards.

How?

- Dispose of unnecessary materials on tables, in drawers, in cabinets and on the floor.
- Remove unnecessary instructions and posters on the walls and on notice boards.
- Attach Red Tags for undecided items for disposal and determine whether the item “may be necessary” or is “unnecessary.”

Hints for actions

- Eliminate all unnecessary items from workplace and reduce clutter.
- Go through all tools, materials, and so forth in the work area.
- Keep only essential items and eliminate what is not required, prioritizing things as per requirements and keeping them in easily-accessible places.
- Everything else should be stored or discarded.

		
<p>Dispose off unnecessary materials on tables, in drawers, in cabinets and on the floor.</p>	<p>The walls and notice boards use X-Y line.</p>	<p>Attach “Red Tags” for undecided items for disposal and determine whether the item “may be necessary” or is “unnecessary.”</p>
		
<p>Good example of chaos free environment at the office desk.</p>	<p>The walls and notice boards use X-Y line.</p>	<p>Remove unnecessary instructions and posters on the wall and on the notice boards.</p>

2S: SETTING/ Material Handling

Basic points for setting material handling

1. Zone/area, clear transport routes;
2. Multi-level shelves and containers;
3. Home for each tool;

How?

- Secure transport routes without obstacles that are even and not slippery.
- Provide multi-level shelves or storage racks near the work area for medical charts and drugs.
- Provide a “home” for medical equipment and work items and mark their position clearly.
- Use carts, hand-trucks and other wheeled devices when moving materials.

Hints for actions

- Organize everything needed in proper order for ease of operation.
- There should be a place for everything and everything should be in its place.
- The place for each item should be clearly labelled or demarcated.
- Items should be arranged in a manner that promotes efficient work flow, with equipment used most often being the most easily accessible.



Clear transport routes with clear mark.



Provide multi-level shelves or storage racks near the work area for medical drugs.



Provide a “home” for each piece of medical equipment.



Use carts, hand-trucks and other wheeled devices when moving materials.



Secure transport routes without obstacles that are even and not slippery.



Provide multi-level shelves for medical charts in the record unit.

2S: SETTING/ Workstation changes

Basic points for setting workstation changes

1. Easy reach

3. Comfortable chairs

2. Elbow height

How?

- Place frequently used tools, controls and materials within easy reach of workers.
- Adjust the working height for each health care worker at elbow level or slightly below it.
- Allow workers to alternate standing and sitting as much as possible and provide good adjustable chairs with good backrests.

Hints for actions


- Items should be arranged in a manner that promotes efficient work flow, with equipment used most often being the most easily accessible.
- Workers should not have to bend repeatedly to access materials.
- Each tool, part, supply, or piece of equipment should be kept close to where it will be used – in other words, straightening the flow path.



Place frequently used tools and materials within easy reach of workers



Adjust the working height for each health care worker at elbow level or slightly below it.



Provide good adjustable chairs with good backrests.



Place frequently used tools within easy reach



Adjust the working height for each health care worker at elbow level or slightly below it to reach.



Set necessary equipment nearby worksite for easy to reach.

2S: SETTING/ Labels, signs, color coding

Basic points for setting labels, signs, color coding


1. Labelling
2. Safety signs
3. Color coding/numbering
4. Signboards/mapping

How?

- Make labels, colors and signs easy to see, easy to read and easy to understand.
- Provide clear instructions and information necessary for safe, efficient work.
- Use color coding rules such as for garbage, medical waste and linen by type.
- Use a board with X-Y axis alignment for sharing medical information, a 5-S corner display and a map of work areas for patients.

Hints for actions

- The place for each item should be clearly labelled or demarcated.
- Items should be arranged in a manner that promotes efficient work flow, with equipment used most often being the most easily accessible.

		
<p>Labels should be easy to read and easy to understand.</p>	<p>Clear instructions for operation of the fire extinguisher.</p>	<p>Numbering rules of medical charts for proper setting.</p>
		
<p>5S corner display by X-Y axis alignment.</p>	<p>The place for each item should be clearly labeled or demarcated.</p>	<p>A map of work areas for patients.</p>

3S: SHINING/ Regular maintenance and cleanliness

Basic points for Shining

1. Cleanliness
2. Waste bin, cleaning tool
3. Proper maintenance of equipment

How?

- Designate the responsibility for keeping areas clean to staff members and clean regularly.
- Place waste bins and cleaning tools in the appropriate workplace.
- Inspect and maintain medical Instruments, machines and tools regularly.

Hints for actions

- Discuss the placement of waste bins in areas where they are needed at your workplace.
- Regular maintenance keeps your medical equipment productive and safe. Check all machine parts carefully. In particular, special care is needed for rotating parts, detachable guards and electrical wires, etc.
- If you find an inexpensive machine, you should check its safety aspects even more carefully. When accidents occur, costs can be enormous.
- Develop a long-term maintenance plan. For instance, you may check the machine thoroughly every weekend.



Designate responsibility for keeping the area clean to staff members.



Place waste bins and cleaning tools in the appropriate workplace location.




Inspect and maintain medical instruments, machines and tools regularly.



Cleaning regularly supports daily work tasks.



Place waste bins at every point where HWs undertake procedures.



Discuss establishing a standard of cleanliness with staff members.

4S: STANDARIZE/ Improving quality of care

Basic points for standardizing



1. Standardized care procedures
2. Procedures of the equipment
3. Ensuring patients' privacy
4. Establishing hand hygiene
5. Reducing patients' waiting time

How?

- Establish the protocols of standardized care procedures about each cure and care, such as maternity, vaccination, medical health checkup.
- Display procedures of the equipment and label clearly what operation is meant.
- Use partitions, curtains and other arrangement for protection privacy of persons cared
- Establish hand hygiene procedures and hygienic washing facilities.
- Reducing waiting time for patients at an outpatient room, a laboratory and a pharmacy.

Hints for actions

- Listing basic clinical procedures makes it easy to find necessary standardized care procedures for documentation.
- Labels and signs are needed for clear instruction for procedures.
- Encouraging collection of patients' opinion is essential for improving quality of care.

 <p>STEPS FOR MALARIA DIAGNOSIS</p> <pre>graph TD A[Malaria Suspected] --> B[Malaria Smear Positive] A --> C[Malaria Smear Negative] B --> D[Non Falciparm] B --> E[Falciparm Species] D --> F[Yes] D --> G[NO] E --> H[.....] C --> I[.....]</pre>	 <p>Instructions for procedures for the use of equipment with operations clearly labeled.</p>	 <p>Curtains for protection privacy of persons cared.</p>
<p>The protocols of standardized malaria diagnosis as a flowchart</p>	<p>Instructions for procedures for the use of equipment with operations clearly labeled.</p>	<p>Curtains for protection privacy of persons cared.</p>
 <p>Proper hand hygiene procedures and hygienic washing facilities.</p>	 <p>A suggestion box helps improve the quality of care.</p>	 <p>Reducing waiting times for patients at outpatient departments, laboratories and pharmacies.</p>
<p>Proper hand hygiene procedures and hygienic washing facilities.</p>	<p>A suggestion box helps improve the quality of care.</p>	<p>Reducing waiting times for patients at outpatient departments, laboratories and pharmacies.</p>

5S: SUSTAIN/ Teamwork

Basic points for sustaining

1. Tool box meeting
2. Record of accidents and injuries
3. Work Improvement Team (WIT), 5S Training

How?

- Hold a brief meeting before beginning work to plan assignments and avoid giving excessive workloads to particular staff members.
- Keep records of medical accidents, occupational injuries and diseases for ensuring both patients and workers' safety.
- Encourage Work Improvement Team (WIT) activities and provide adequate participatory 5S and safety/health improvement training to all workers.

Hints for actions

- Record all accidents, absences, sickness, and other health-related events.
- Plan annual work schedules including sufficient training periods.



Hold a brief meeting before work



Keep records of both medical accidents and occupational injuries



Visit 5S team regularly at the workplace



5S training for all staff members



Holding monthly WIT meetings helps maintain 5S activities



5S facilitators organize engaging 5S workshops.

ACTION PLAN FOR CQI/KAIZEN ACTIVITY

CQI topic: **Making bigger workshop space for start maintaining electronic equipment in dust-free environment**

Date: 22 / 10 / 2018

Place: **** Workshop

Implementer(s): Workshop staff

Annex 6(b): CQI Action Plan Form (SAMPLE)

No.	Actions to be taken	From: 22/10/2018 to: 21/11/2018			Responsible person(s)	Resources	Tick when completed
		1M	2M	3M			
1	Get permission to open the two un-used offices in the Workshop building	●			Manager		✓
2	Open the locks of two offices	●			Denis	Tools for break door locks	✓
3	Clean the two offices	→			All cleaners		✓
4	Discuss and agree on needs and layout	●			All staff		✓
5	Relocate desk and chairs	●			Moris		✓
6	Acquire 1-2 shelves from the hospital	→			Manager & all staff	Shelves	✓
7	Transfer tools and cables	●			Daniel		✓
8	Transfer electronics equipment	●			Bonny		✓
9	Start maintenance of electronics equipment in dust-free room	↑			All staff		✓
10							

COMPLETED

After filling in the above, everyone signs an approval.

Ogwai Walter	Akana Moris	S. Mulepo
Obia Daniel	Doreen Mubiru	Naoki Mimuro

Annex 7: Quarterly Report Format

(RRH/RW Headed Paper)

Date: (Date/Month/Year)

To: CHS (Health Infrastructure)
Health Infrastructure Department, MOH

Thru: Director/RRH

Re: "X" Regional Medical Equipment Maintenance Workshop Quarter "Y"
Performance Report

1. Executive Summary

This is a summary report of maintenance activities carried out by "Name of Workshop" for the period **October – December 2019**. The report highlights the progress of implementation of planned activities, budget performance, cost of work undertaken, cost of spare parts for pending work, Medical equipment inventory operational status, challenges, Work plan for the next quarter, concluding remarks and recommendations.

For the period under review, the workshop performance improved and the following achievements were registered:

- 1) Medical equipment in good functional condition increased from 69% to 73%.
- 2) To date medical equipment inventory update in NOMAD stands at 100% (RRH), 95% (GHs), 80%(HCIVs) and 35% (HCIIIs).
- 3) 12 Midwives were trained in operation, care and basic maintenance of Oxygen concentrators, Suction pump and Baby incubators from 8HCIIIs.
- 4) The Budget performance was 87%. The balance of the funds is committed but await delivery of office stationery and supplies by the service provider in 1 to 2 weeks' time.
- 5) 50 oxygen cylinders were produced and delivered to Jinja RRH (24 cylinders), 2 GHs (14 cylinders) and 4 HCIVs (12 cylinders).

The main challenges were:

- 1) Delayed release of funds to the RW team and this resulted in delayed implementation of planned activities.
- 2) The procurement of spare parts was very slow and this resulted in long downtime for some of the medical equipment.

To overcome the above challenges, the management of the RRH should facilitate fast release of funds for planned activities and secure framework contracts for supply of fast-moving spare parts.

The detailed performance report for quarter “Y” is outlined below.

2. Planned activities and progress of implementation

Planned maintenance activities were carried out in the following health facilities during the period under review.

RRH	Mazuungu
GH	Maizimasa, Kibaale, Kiboko
HC IV	Ma-aji, Ndeku, Mukokolo, Alia

The status of implementation of the planned activities is summarized in the Table below:

S/N	Planned Activity	Progress of Implementation	Achievements/Outputs & Immediate outcomes	Actual expenditure for Activities (UGX)	Remark
	Carry out medical equipment servicing and repair in 11 HFs	Completed	315 pieces of equip. were maintained and 11 are pending	3,450,000	Main focus was servicing Lab equipment
	Carry out Medical Equipment inventory collection and update for ** GH and ** HCIV	Ongoing	567 pieces of equipment entered in NOMAD	1,450,000	Activity was funded by IP
	Conduct user training on operation & maintenance of concentrators and incubators for 16 nurses from 8 HCIIIs	Completed	12 Midwives trained	2,000,000	2 of the target HCIIIs had only one midwife instead of 2
	Attend the performance review meeting for RWs	Completed	Successful meeting held in Pinga RRH	1,300,000	HD & 2 BMETs attended,
	5S-CQI-TQM – Spare parts management in Workshop Store	Completed	Procured stock cards, store issue forms and shelves	260,000	New shelves installed & store space increased

S/N	Planned Activity	Progress of Implementation	Achievements/Outputs & Immediate outcomes	Actual expenditure for Activities (UGX)	Remark
	Oxygen Plant operation and production	Ongoing	50 cylinders produced and delivered to RRH, 2 GHs and 4 HCIV.	-	Service contract will be expired in Sep. 2019.

Detailed description of the equipment maintenance activities is in the attached Work Record sheets – *Annexes 7(a) & 7(b)*.

3. Financial expenditure and budget performance for the quarter

The table below summarizes the expenditure for the period as extracted from Annexes 1 and 2.

Reference	Item Description	Amount (UGX)
Quarterly Release	Allocation for the quarter	
-	Balance brought forward from previous quarter	
-	Additional fund from IP(s)	
Total Available Budget for Quarter		
-Spare parts	Cost of spare parts purchased in the quarter	
-Per diem & SDA	Per diem & SDA expenditure for the quarter	
Fuel/Mileage Expenditure	Mileage/Fuel: Km (Litres of Diesel)	
Other Workshop Operational costs		
1) Electricity bill	Electricity cost	
2) Water bill	Water cost	
3) Office stationery	Cost of office stationery purchased in the quarter	
4) Vehicle maintenance	Cost of vehicle maintenance/repair	
5) Technical Supervision and Coordination	RW Management Committee Meeting/RW performance review meeting cost	
6) Capacity building training	Cost of User training for Nurses and Techs.	
Total Expenditure (UGX)		
Overall Budget Performance (Total Expenditure/Total Available Budget for Quarter) %		

4. Medical equipment inventory analysis report for the quarter

Based on the updated medical equipment condition in the NOMAD database, the equipment condition was as shown in the table below.

This Quarter

Equipment condition	A	B	C	D	E	F	Total
Quantity							
Percentage (xx.x%)	%	%	%	%	%	%	100.0%

Previous Quarter

Equipment condition	A	B	C	D	E	F	Total
Quantity							
Percentage (xx.x%)	%	%	%	%	%	%	100.0%

This analysis is based on the medical equipment inventory for 1RRH, 5GHs and 11HCIVs.

5. Challenges

During implementation of the planned activities, the following challenges were encountered:

- 1) There was delayed implementation of planned activities due to late release of funds.
- 2) The procurement of spare parts delayed and this resulted in long downtime for some of the medical equipment.

To overcome the above challenges, the management of the RRH should facilitate fast release of funds for planned activities and secure framework contracts for supply of fast-moving spare parts.

6. Proposed work plan and budget for the next quarter

The work plan and budget for the next quarter is attached for perusal [*Annex 7(d)*]. The list and estimated cost of required spare parts for pending work is attached as *Annex 7(c)*. Any emergency repair shall be handled on a case by case basis. The implementation schedule for the planned maintenance activities is attached as *Annex 7(e)*.

7. Recommendations and Concluding Remarks

There was marked improvement in the equipment functional condition as a result of repair of a number of equipment that was pending for the last 3 months.

In order to reduce on the number of equipment that require repair, the Top Management will need to prioritize procurement of spare parts and training of the equipment users.

Prepared and signed by: _____
"Name and Signature of Workshop Manager"

C.c DHOs of the Districts in the Catchment Area

C.c All Medical Superintendents of Health Facilities in Catchment Area (i.e. GHs and HCIVs)

ANNEX 7(a): WORK RECORD SHEET

Item No.	Health Facility Name	Number of job cards raised	Number of Days spent at Health Facility	Total cost of spare parts used	Total number of pending jobs	Estimated cost for spare parts for pending jobs	Number of BME/Ts involved in maintenance activity
1	*** GH	30	2 days	500,000	3	300,000	3
2							
3							
4							
Total							

Prepared by: _____ Date: _____
 "Name and Signature of Workshop Manager"

Reviewed and approved by: _____ Date: _____
 "Name and Signature of HD/ RRH"

ANNEX 7(b): HEALTH FACILITY JOB CARD RECORD SHEET

Job card serial No.	Health Facility Name	Equipment Name	Dept/ Section	Description of spare part used	Qty	Final Job Status (<i>Complete/Pending</i>)	Spare parts required for completing pending work	Date of next planned maintenance visit
025	*** HCIV	Operation table	OPD	Hydraulic oil	2Lts	Complete	NIL	Aug. 2019
026	**** GH	Operation light	Surgey	Surgical Lamp, 20W LED	02	Complete	NIL	Oct. 2019
027	*** GH	Microscope	Lab	-	-	Pending	Bulb, 9V	Oct. 2019

ANNEX 7(d): PROPOSED WORK PLAN AND BUDGET FOR QUARTER “X”

S/ N	Planned Output Target	Planned Activities to Deliver Outputs	GOU Quarterly allocation	IP Quarterly allocation	TOTAL
			Q3	Q3	
1	Carry out PPM and PCM in 2RRHs, 8GHs, 46HCIVs and at least 60HCIIIs in ** region	PPM and PCM in hospitals, HCIVs and HCIIIs in ** region	80,000,000	20,000,000	100,000,000
2	Procure assorted Medical Equipment spare parts worth UGX 285m	Procurement of assorted Medical Equipment spare parts	53,000,000	-	53,000,000
3	Train 25 BME/Ts	Biomedical Engineering training in maintenance of ICU, Oxygen therapy and theatre equipment	40,000,000	30,000,000	70,000,000
4	Collect and update Medical equipment in 90% of the HFs in catchment area	Collection and update of Medical equipment inventory	15,000,000	20,000,000	35,000,000
5	Attend 4 quarterly Regional Workshop performance review meetings	Attend regional equipment Workshops' performance review meetings	1,000,000	4,000,000	5,000,000
6	Calibrate all the tool kits for servicing auxiliary lab equipment	Carryout calibration of tool kits for servicing auxiliary laboratory equipment	-	60,000,000	120,000,000
7	Emergency work		3,000,000	-	3,000,000
TOTAL			188,000,000	134,000,000	322,000,000

Prepared by: _____ Date: _____
 “Name and Signature of Workshop Manager”

Reviewed and Approved by: _____ Date: _____
 “Name and Signature of HD/RRH”

Annex 8: Job Card Form

JOB CARD

Date: / /

Job Card Serial Number:

Health Facility	Department	Section/ Room
Name of Equipment/ Furniture/ Other repair works		Original Operational Status A. B. C. D. E. F.
Serial No./ Engraved No.		Model Name / Model No.
Fault Reported/ Detected		
Work done		
Extra work to be done/Advise to the Health Facility/Remarks		
Name of Eng/Tech		Signature
Name of Casual Labourer	Contact No.	Cost
Spare Parts Used and Other Supplies		
Name		Qty
Final Operational Status After Maintenance: A. B. C. D. E. F. (tick or circle)		
Estimate Date of Next Service (Month/Year)		/
Follow-up Contact No. for Maintenance Team Member		Tel
Comment and Endorsement by User Department/Health Facility in charge		
Endorsed by	Name	Stamp
	Title	
	Signature	
	Date	

Note; Original- kept by Workshop, 1st Copy- kept by HF, 2nd Copy-Retained in Job Card Book

Annex 9: Annual Assessment Sheet for RWs

Annual Assessment for the RWs, Year _____
 Scoring Criteria: 3: Good, 2: Fair, 1: Poor/low, 0: Not Performed

No	PERFORMANCE INDICATORS	Central	Mbale	Soroti	Lira	Gulu	Arua	Hoima	Fort Portal	Kabale	Mubende	Moroto	Jinja	Masaka	Mbarara	Naguru	Average score
1	Available WS staff for ME maintenance																
2	Timely release of WS budget to WS team and allocative efficiency																
3	ME inventory properly updated and data entered in NOMAD * <i>double score</i> (full mark: 6)																
4	Job cards properly prepared and used																
5	Productivity: No. of job cards raised / No. staff involved per quarter >=50																
6	Adequate spare parts purchased in timely manner																
7	Planned Preventive Maintenance periodically carried out for Lab. equipment and other selected equip																
8	Routine maintenance carried out at least once a quarter for all hospitals and HCIs																
9	User training planned for and conducted using WS budget																
10	Availability of integrated Workplan and budget developed by WS staff and RRH management																
11	Timely submission and quality of Quarterly WS progress Reports																
12	Teamwork/ Team building																
13	Continuous implementation of CQI activities in the WS																
		TOTAL (Full-mark:42)															
		TOTAL % (score/42 x 100)															
		Previous TOTAL %															
		Increment (%) from previous score															

Annex 10: Maintenance Requisition Form

***** Regional Medical Equipment Maintenance Workshop

MAINTENANCE REQUISITION FORM

Requisition Form Serial No:

Health Facility	Department	Section/ Room
Name of Equipment/ Furniture/ Others	Model Name / Model No.	Equipment Serial No./ Engraved No.

Description of the problem identified

Details of Person Reporting

Name:	Title:	Date:
-------	--------	-------

Endorsed by Health Facility in charge

Received by member of the Workshop

Name	Name
Title	Title
Signature & Stamp	Signature
Date	Contact No. Date

Note; Original- kept by Workshop, Copy- kept by Health Facility/Health Unit

Annex 11: Maintenance Service Sticker Form

(80 x 50)mm

 **FORT-PORTAL REGIONAL MEDICAL
EQUIPMENT MAINTENANCE WORKSHOP**

JOB CARD NO.

DATE OF SERVICES:

NEXT SERVICE: 3M 6M 12M

CURRENT CONDITION: A B C D E F

SAFETY TEST USER TRAINING

SERVICE CONTRACT: YES NO

NAME SIGN.....

CONTACT:



Supported by

